

NDI Passive Spheres™


Single-use NDI Passive Spheres™ are used on optically-tracked tools during computer-assisted therapies including image-guided surgery.

Packaged for ease-of-use in a sterile environment, NDI Passive Spheres™ are designed by the same company that provides the optical measurement technology at the core of computer-assisted therapy systems used in operating rooms around the world since 1997.



NDI Passive Spheres™ offer ...

- Single-use, sterile packages of 1, 3 or 4 units to eliminate waste, costs and degradation associated with re-processing
- A reliable snap-fit to NDI mounting posts for a secure attachment to tools
- Optimal design by NDI for accurate tracking
- Convenient 11.5 mm (0.45 ") diameter for best results and minimal obstruction
- Exclusive distribution through your Scanlan International distributor representative to assure local stock, service and support.



Established more than 20 years ago, NDI provides the core technology and related accessories for advanced 3D measurement and tracking applications used by many of the world's leading hospitals and medical research institutions.

In computer-assisted therapies from neurosurgery to orthopedic surgery, NDI technology is trusted for its accuracy and reliability. Today, the company supports over 5,000 installations in more than 25 countries around the world.

IMAGE GUIDED SURGERY

- 8800737 (1) Marker per pack, (5) Packs per box
- 8800738 (3) Markers per pack, (5) Packs per box
- 8800739 (4) Markers per pack, (5) Packs per box

NDI is a registered trademark of Northern Digital Inc.



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

ONE SCANLAN PLAZA • SAINT PAUL, MN 55107 • USA

800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com

© 2005 Scanlan International, Inc.