

VATS MICS Wolf Suction Instrument – 3.2 mm
Developed in cooperation with Randall K. Wolf, MD

Working lengths:
Pediatric/Small Adult 6-3/4" / 17 cm
Adult 8-3/4" / 22 cm
Shaft Diameter: 3.2 mm
Shaft: Curved
Variable Suction Option: Yes
Tip Width: 6 mm



Enlarged To Show Detail

VATS MICS Dennis Suction Instrument – 5 mm

Working lengths:
Adult 8-3/4" / 22 cm
Adult 12-1/2" / 32 cm
Shaft Diameter: 5 mm
Shaft: Curved
Variable Suction Option: Yes
Tip Width: 7.5 mm



Enlarged To Show Detail

VATS D'Amico Suction Instrument – 5 mm
Developed in cooperation with Thomas A. D'Amico, MD

Working lengths:
Adult 12-1/2" / 32 cm
Shaft Diameter: 5 mm
Shaft: Straight
Variable Suction Option: Yes
Tip Width: 5 mm



Enlarged To Show Detail

VATS MICS D'Amico Suction Instrument – 10 mm
Developed in cooperation with Thomas A. D'Amico, MD

Working lengths:
Adult 12-1/2" / 32 cm
Shaft Diameter: 10 mm
Shaft: Straight
Variable Suction Option: Yes
Tip Width: 10 mm




Enlarged To Show Detail

VATS / MICS Suction Instruments

Catalog Number	Description	Shaft Diameter	Working Length	Total Length
Wolf Suction Instruments				
9009-908	Curved Blunt Tip, 6 mm width	3.2 mm	6-3/4" / 17 cm	12-1/4" / 31 cm
9009-909	Variable Suction, Curved Blunt Tip, 6 mm width	3.2 mm	6-3/4" / 17 cm	12-1/4" / 31 cm
9009-910	Curved Blunt Tip, 6 mm width	3.2 mm	8-3/4" / 22 cm	14-1/4" / 36 cm
9009-911	Variable Suction, Curved Blunt Tip, 6 mm width	3.2 mm	8-3/4" / 22 cm	14-1/4" / 36 cm
Dennis Suction Instruments				
9009-914	Curved Blunt Tip, 7.5 mm width	5 mm	8-3/4" / 22 cm	14-1/4" / 36 cm
9009-915	Variable Suction, Curved Blunt Tip, 7.5 mm width	5 mm	8-3/4" / 22 cm	14-1/4" / 36 cm
9009-916	Curved Blunt Tip, 7.5 mm width	5 mm	12-1/2" / 32 cm	17-3/4" / 45.5 cm
9009-917	Variable Suction, Curved Blunt Tip, 7.5 mm width	5 mm	12-1/2" / 32 cm	17-3/4" / 45.5 cm
D'Amico Suction Instruments				
9009-900	Straight blunt tip	5 mm	12-1/2" / 32 cm	17-3/4" / 45.5 cm
9009-901	Variable suction, straight blunt tip	5 mm	12-1/2" / 32 cm	17-3/4" / 45.5 cm
9009-920	Straight blunt tip (for decortication)	10 mm	12-1/2" / 32 cm	17-3/4" / 45.5 cm
9009-921	Variable suction, straight blunt tip (for decortication)	10 mm	12-1/2" / 32 cm	17-3/4" / 45.5 cm
Suction Accessories				
9009-905	Brush set, pack of 10, for D'Amico, Wolf, Dennis suction instruments			16-1/2" / 42 cm
9009-922	Brush set, pack of 10, for 10 mm D'Amico suction instruments			19" / 48.3 cm

For more information or to place an order,
please contact your local Scanlan representative, or call...
U.S. & Canada: 800-328-9458
International: 651-298-0997
FAX: 651-298-0018
www.scanlaninternational.com




Surgical products designed and manufactured by the Scanlan family since 1921
ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com

MK-B-VMSI Rev. A
© 2011 Scanlan International, Inc.

VATS / MICS Suction Instruments



Surgical products designed and manufactured by the Scanlan family since 1921



VATS MICS Suction Instruments



VATS MICS Wolf Suction Instrument – 3.2 mm
Developed in cooperation with Randall K. Wolf, MD

- Tip is blunted so it can safely be used as a dissector
- Option of 17 cm or 22 cm working lengths for pediatric / adult patients
- Yankauer handle provides familiar, secure grip




9009-908
Wolf Suction Instrument (Actual Size)

9009-910
Wolf Suction Instrument (Actual Size)

VATS MICS Dennis Suction Instrument – 5 mm

- 5 mm shaft and tip opening won't clog during procedure
- Tip is blunted so it can safely be used as a dissector
- Option of 22 cm or 32 cm working lengths
- Option of with or without fingertip controlled suction regulator



9009-914
Dennis Suction Instrument (Actual Size)

9009-916
Dennis Suction Instrument (Actual Size)

VATS D'Amico Suction Instrument – 5 mm
Developed in cooperation with Thomas A. D'Amico, MD

- Specifically developed for VATS procedures for both suction and dissection
- Tip is blunted so it can be safely used as a dissector
- Option of with or without fingertip controlled suction regulator



9009-900
D'Amico Suction Instrument (Actual Size)

VATS MICS D'Amico Suction Instrument – 10 mm
Developed in cooperation with Thomas A. D'Amico, MD

- VATS drainage of empyema and decortication procedures
- Drainage of large hematoma and clots
- Option of with or without fingertip controlled suction regulator



9009-920
D'Amico Suction Instrument (Actual Size)