

# EASY-TAG™ Tracking System

New PINK & GREEN Colors

- 
- Carts, shelves, and more
  - Endoscopes
  - Instrument trays
  - Instrument repairs
  - Trays with missing instruments
  - Vendor / loaner trays

SCANLAN ★

INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*

# EASY-TAG™ Tracking System

*Designed to answer the need for a cost effective, easy, secure method of tracking instruments, trays, equipment, and repairs through the sterilization, decontamination wash and cleaning process.*


Tracking Repairs



Vendor Loaner Trays



Tracking Endoscopes



Missing Instruments



- **Sterilizable** - EO, Steam, and Gas Plasma
- **Versatile Surface** - allows placement of pre-printed labels or can be written on both sides of tag
- **No Adhesives** - quick, easy and secure
- **Easy to Apply** - the flexible handle loops around instruments and tray handles
- **Easy to Remove** - Tear-Away™ handle provides easy separation after use
- **Color Coded** - for fast and easy identification
- **Latex-free**



<b>1001-904</b>	EASY-TAG, Black
<b>1001-905</b>	EASY-TAG, Red
<b>1001-906</b>	EASY-TAG, Yellow
<b>1001-907</b>	EASY-TAG, Blue
<b>1001-908</b>	EASY-TAG, Green
<b>1001-909</b>	EASY-TAG, Pink

250 Tags per spool  
2.5" x 4.75" / 63.5 mm x 120.65 mm  
Elastomer Handle Polymer Tag

**For more information or to place an order,**  
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

[www.scanlaninternational.com](http://www.scanlaninternational.com)

**SCANLAN** ★

INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.  
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • [www.scanlaninternational.com](http://www.scanlaninternational.com)

MK-B-EZPG  
© 2010 Scanlan International, Inc.