

# Specialty PEDIATRIC INSTRUMENTS



INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*



Specialty

# PEDIATRIC INSTRUMENTS

Scissors, Clamps, Needle Holders, Retractors, Forceps

Since 1921, Scanlan International, Inc. has been designing and manufacturing the highest quality surgical instrumentation. Our pediatric line is a unique collection of instruments for procedures with small to medium sized surgical fields. Each instrument is designed with the highest grade stainless steel or titanium alloys, and are handcrafted with meticulous precision... allowing excellent tactile response and control.

# PEDIATRIC INSTRUMENTS

Specialty



Scissors, Clamps, Needle Holders, Retractors, Forceps

This collection has been selected from our full line of high quality surgical instrumentation (consisting of over 3000 designs). Please visit us at [www.scanlaninternational.com](http://www.scanlaninternational.com) to review additional selections and options. Or simply scan the QR code with your smart phone to connect you directly to our site.



## NEEDLE HOLDERS



6006-58

NEEDLE HOLDERS Stainless Steel	Cat. No.	Description	Length	Suture Size**
	6006-44	Round, spring, with ratchet	4 3/4" (12 cm)	7-0
	6006-58	Flat, spring, with ratchet	5 1/2" (14 cm)	5-0
	6006-50	Round, spring, with ratchet	5 3/4" (14.5 cm)	5-0
	6006-29	Ring, Ryder-style, micro	6 1/4" (16 cm)	5-0
	6006-32	Ring, micro	7 1/4" (18.5 cm)	7-0
	6006-120	Round, micro, with ratchet	8 1/4" (21 cm)	5-0
	6006-36	Ring, micro	8 3/4" (22 cm)	5-0
	6006-40	Ring, micro	10" (25.4 cm)	5-0



3003-430

NEEDLE HOLDERS Titanium	Cat. No.	Description	Length	Suture Size**
	3003-90A	Round, micro, with ratchet	4 3/4" (12 cm)	5-0
	3003-207	Flat, spring	5 1/2" (14 cm)	7-0
	3003-430	Ring, Ryder-style	6" (15 cm)	5-0
	* 3003-450	Ring, micro	6" (15 cm)	6-0
	3003-93A	Round, micro, with ratchet	6" (15 cm)	5-0
	3003-432	Ring, Ryder-style	6 3/4" (17 cm)	5-0
	3003-530	Round, with ratchet	9" (23 cm)	2-0

\*\*Suture size listed or smaller

## FORCEPS



4004-56

FORCEPS Stainless Steel	Cat. No.	Description	Length
	4004-40	Flat, DeBakey, 1 mm tips	6" (15 cm)
	4004-42	Flat, DeBakey, 2 mm tips	6" (15 cm)
	4004-56	Flat, gold carbon inlay tips	7" (18 cm)
	4004-230	Round, DeBakey	7 1/8" (18.1 cm)
	5005-33	Ring, curved jaws	7 1/4" (18.5 cm)



3003-150

FORCEPS Titanium	Cat. No.	Description	Length
	3003-150	Round, platform jaw	4 3/4" (12 cm)
	3003-154	Round, micro	6" (15 cm)
	3003-03	Flat, DeBakey	6" (15 cm)
	3003-730	Round, DeBakey	7 1/8" (18.1 cm)
	3003-05	Flat, DeBakey	7 3/4" (19.5 cm)
	3003-778	Round, 1 x 2 micro teeth	8 1/4" (21 cm)

## SCISSORS



7007-420

SCISSORS Spring, Round Handle	Cat. No.	Description	Length
	7007-420	25° angled blades	4 1/2" (11.5 cm)
	7007-422	45° angled blades	4 1/2" (11.5 cm)
	7007-108	Curved blades, sharp tips	5 1/2" (14 cm)
	7007-590	25° angled blades	6 3/4" (17 cm)
	7007-592	45° angled blades	6 3/4" (17 cm)
	7007-594	60° angled blades	6 3/4" (17 cm)



7007-410

SCISSORS Spring, Flat Handle	Cat. No.	Description	Length
	7007-401N	25° angled blades	4 1/8" (10.5 cm)
	7007-410	45° angled blades	4 1/2" (11.5 cm)
	7007-412	Micro blades, sharp tips	4 1/2" (11.5 cm)
	7007-539N	25° angled blades	6" (15 cm)
	7007-452	45° angled blades	6 1/2" (16.5 cm)

## SCISSORS



7007-210SC

SCISSORS Ring Handle	Cat. No.	Description	Length
	7007-210SC	Tenotomy tips	4 5/8" (11.8 cm)
	7007-214SC	Curved blades	4 3/4" (12 cm)
	* 7007-395SC	Sternum Pediatric Super Cut™ Scissors	5" (12.7 cm)
	7007-216SC	Curved blades	6" (15 cm)
	7007-263SC	Tenotomy tips	6 1/2" (16.5 cm)
	7007-262SC	SCANLAN® Puskas™ "Black Knight", Sharp Tips	6 1/2" (16.5 cm)
	7007-264SC	Tenotomy tips	7" (18 cm)
	7007-266SC	Sharp tips	7" (18 cm)
	7007-211SC	Curved blades	7" (18 cm)
	* 7007-397SC	Sternum Pediatric Super Cut™ Scissors	8" (20 cm)
	7007-222-2SC	Curved blades	9" (23 cm)
	7007-224SC	Curved blades	10 1/4" (26 cm)

\*Custom made at time of order. Firm sale. Custom modifications available upon request.



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921



Specialty

# PEDIATRIC INSTRUMENTS

Scissors, Clamps, Needle Holders, Retractors, Forceps

## CLAMPS



2002-235

CLAMPS  
Stainless Steel

Cat. No.	Description	Length
2002-236	60° angled jaws	5 1/4" (13.5 cm)
2002-234	25° angled jaws	5 3/4" (14.5 cm)
2002-235	45° angled jaws	5 3/4" (14.5 cm)
2002-229	Small spoon jaws, straight shanks	5 3/4" (14.5 cm)
5005-60	Bailey, 30° angled jaws	7" (18 cm)
5005-62	Bailey, 60° angled jaws	7" (18 cm)
5005-64	Bailey, 90° angled jaws	7" (18 cm)
* 2002-51	DeBakey, straight jaws	12 1/2" (32 cm)



3003-620

CLAMPS  
Titanium

Cat. No.	Description	Length
3003-616	6 mm jaws, slightly angled shanks, 14 mm jaw length	4 1/8" (10.5 cm)
3003-618	6 mm spoon jaws, straight shanks, 14 mm jaw length	4 1/8" (10.5 cm)
3003-620	9 mm jaws, slightly angled shanks, 16 mm jaw length	4 1/8" (10.5 cm)
3003-614	60° angled jaws, 23 mm jaw length	4 1/8" (10.5 cm)
3003-613	45° angled jaws, 25 mm jaw length	4 1/4" (11 cm)
3003-634	60° angled jaws, 30 mm jaw length	4 1/4" (11 cm)
3003-628	Left angled jaws, 23 mm jaw length	4 1/4" (11 cm)
3003-630	Right angled jaws, 23 mm jaw length	4 1/4" (11 cm)
3003-633	45° angled jaws, 30 mm jaw length	4 1/4" (11 cm)
* 3003-636	6 mm jaws, slightly angled shanks, 22 mm jaw length	4 3/8" (11.1 cm)
3003-637	9 mm jaws, slightly angled shanks, 20 mm jaw length	4 3/8" (11.1 cm)
3003-72	Lambert-Kay aorta clamp	8" (20 cm)

## RETRACTORS



1001-10

RETRACTORS  
Hand-Held

Cat. No.	Description	Length
1001-10	3/4" wide	6 1/2" (16.5 cm)
1001-12	1/2" wide	6 1/2" (16.5 cm)
1001-14	1/4" wide	6 1/2" (16.5 cm)



8008-43

RETRACTORS  
Sternal

Cat. No.	Description
* 8008-41	Finochietto Neonate Sternal Retractor, aluminum, 36 mm spread; includes two blades: each 7 mm long x 7 mm deep; two arms: each 37 mm long
* 8008-42	Finochietto Neonate Sternal Retractor, titanium, 36 mm spread; includes two blades: each 7 mm long x 7 mm deep; two arms: each 37 mm long
* 8008-43	Titanium Retractor with interchangeable blades, 52 mm spread; includes two pair of blades: one pair 10 mm long x 15 mm deep, one pair 20 mm long x 15 mm deep; two arms: each 40 mm long
* 8008-45	Bailey Finochietto Sternal Retractor, titanium, 80 mm spread; includes two blades: each 20 mm long x 20 mm deep; two arms: each 57 mm long
* 8008-47	Bailey Finochietto Retractor, aluminum, 80 mm spread; includes two blades: each 15 mm long x 15 mm deep; two arms: each 56 mm long

For more information or to place an order, please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458  
International: 651-298-0997  
FAX: 651-298-0018  
[www.scanlaninternational.com](http://www.scanlaninternational.com)



\*Custom made at time of order. Firm sale.  
Custom modifications available upon request.



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921  
ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.  
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • [www.scanlaninternational.com](http://www.scanlaninternational.com)

MK-B-PEDINST Rev. A  
© 2012 Scanlan International, Inc.