

SCANLAN®  
Single-Use & Reusable  
**Aorta/Vein Punches**



INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*

# AORTA / VEIN PUNCH

Single-Use


SCANLAN®

## Single-Use Aorta/Vein Punch

- SCANLAN® Scissor Cut™ blades featured on the aorta / vein punches provide clean accurate incisions.
- **Tapered tip** for an unobstructed view of the operating field.
- **Available in two lengths:**  
Original: 5 3/4" (14.5 cm)  
Ultimate: 7 1/2" (19 cm)
- **Available in nine cutting sizes**  
7 Round  
2 Elliptical
- **Radiopaque**
- **Not made with natural rubber latex.**




The SCANLAN® Single-Use Aorta /Vein Punch tip is inserted into a small incision in the vessel wall...



The double action Scissor Cut® blade is drawn upward through the tissue...



...Providing a clean, accurate incision, leaving the anastomotic site virtually undamaged.




Use your Smart Phone and QR Reader software to photograph the code which will demonstrate the use of the Single-Use Aorta Vein Punch.



# ORIGINAL Aorta / Vein Punch

Single-Use



SCANLAN®

Original Aorta / Vein Punch

Cat. No.

Description

Total Length

		Cat. No.	Description	Total Length
Round		1001-600	● 2.7 mm	5-3/4" (14.5 cm)
		1001-601	● 3.5 mm	
		1001-602	● 4.0 mm	
		1001-604	● 4.5 mm	
		1001-605	● 5.0 mm	
		1001-606	● 5.5 mm	
		1001-607	● 6.0 mm	
Elliptical		1001-608	● 2.0 x 4.0 mm	
		1001-610	● 2.0 x 6.0 mm	

1 per sterile package, 5 packages per box





[www.scanlaninternational.com](http://www.scanlaninternational.com) | [facebook](#) | [twitter](#) | [instagram](#)

SCANLAN INTERNATIONAL, Inc. • One Scanlan Plaza • Saint Paul, Minnesota 55107 • U.S.A. • 800-328-9458 • 651-298-0997 • FAX 651-298-0018  
©2016 Scanlan International, Inc.

# ULTIMATE Aorta / Vein Punch

Single-Use



SCANLAN® Ultimate Aorta / Vein Punch		Cat. No.	Description	Total Length
Round		1001-620	● 2.7 mm	7-1/2" (19 cm)
		1001-621	● 3.5 mm	
		1001-622	● 4.0 mm	
		1001-624	● 4.5 mm	
		1001-626	● 5.0 mm	
		1001-628	● 5.5 mm	
		1001-629	● 6.0 mm	
Elliptical		1001-609	● 2.0 x 4.0 mm	
		1001-611	● 2.0 x 6.0 mm	

1 per sterile package, 5 packages per box



# AORTA / VEIN PUNCH


Reusable

SCANLAN®


Reusable Aorta/Vein Punch

Precisely engineered Scissor

Cut™ blades



New Tray  
2135-04



Actual Size

Reusable Aorta / Vein Punch  
Round

SCANLAN®

Reusable Aorta / Vein Punch

Cat. No.

Description

Total  
Length

		Cat. No.	Description	Total Length
Round		1001-03	● 2.7 mm With Sterilizing Box & Tray (2135-04)	5-3/4" (14.5 cm)
		1001-04	● 4.0 mm With Sterilizing Box & Tray (2135-04)	
		1001-05	● 4.5 mm With Sterilizing Box & Tray (2135-04)	
		1001-06	● 5.5 mm With Sterilizing Box & Tray (2135-04)	

NEW SCANLAN® V.I.P.™

Metal Sterilization Tray

Cat. No.

Description

	2135-04	Sterilizing Box & Tray
--	---------	------------------------



[www.scanlaninternational.com](http://www.scanlaninternational.com) | [facebook](#) | [twitter](#) | [instagram](#)

SCANLAN INTERNATIONAL, Inc. • One Scanlan Plaza • Saint Paul, Minnesota 55107 • U.S.A. • 800-328-9458 • 651-298-0997 • FAX 651-298-0018  
©2016 Scanlan International, Inc.

# AORTA / VEIN PUNCH

Reusable

SCANLAN®  
Reusable Aorta/Vein Punch  
Assembly diagram





## Markers Designed For Assessment of Graft Patency During Postoperative Angiography

### SCANLAN® Radiopaque Graft Markers

Scanlan radiopaque graft markers are used to permanently mark the proximal anastomosis during coronary or peripheral bypass graft procedures for postoperative angiography assessment. Made from a lightweight flexible silicone material that minimizes reflective glare during a CT scan or MRI.

### SCANLAN® Radiopaque Graft Markers

	Cat. No.	Description
	1001-83	A/C LOCATOR® graft marker
	1001-831	RADIOMARK® graft marker

Not made with natural rubber latex.  
1 per sterile package, 10 packages per box.

For more information or  
to place an order,  
please contact your local  
Scanlan representative,  
or call...

U.S. & Canada: 800-328-9458  
International: 651-298-0997  
FAX: 651-298-0018  
[www.scanlaninternational.com](http://www.scanlaninternational.com)

