

SCANLAN® RELIANCE Bulldog Clamps



Actual Size




INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

SCANLAN® Reliance Bulldog Clamps

Reliable, Surgeon-Controlled Vessel Clamps

The Scanlan® Reliance Bulldog Clamp has been designed to provide secure, robust clamping pressure to assist in the temporary occlusion of vessels. Manufactured with the highest quality stainless steel, its 1 x 2 DeBakey jaws ensure a firm yet atraumatic clamping pressure. A unique suture port allows optional suture tagging to enhance surgical location and recovery.



Features Include...

- 1 x 2 DeBakey jaws
- 350 grams closing pressure*
- Suture port
- Highest quality stainless steel
- Reusable / durable
- Low profile design
- Radiopaque

Reliance Bulldog Applier / Introducer

- Custom jaw allows for secure deployment of bulldog into the surgical field
- Designed for use with 12 mm cannula

* Approximate clamping pressure



SCANLAN® Reliance Bulldog Clamps

Reliable, Surgeon-Controlled Vessel Clamps

SCANLAN® Reliance Bulldog Clamps Straight

	Description
3795-52	40 mm Straight DeBakey Jaws 2 3/4" (7 cm)
3795-50	25 mm Straight DeBakey Jaws 2" (5 cm)
3795-56 NEW	15 mm Straight DeBakey Jaws 1 5/8" (4 cm)
3795-54	7 mm Straight DeBakey Jaws 1 1/4" (3.2 cm)

SCANLAN® Reliance Bulldog Clamps Curved

	Description
3795-51	25 mm Curved DeBakey Jaws 2" (5 cm)
3795-57 NEW	15 mm Curved DeBakey Jaws 1 5/8" (4 cm)

Firm Sale - May Be Available For Evaluation

SCANLAN® Reliance Bulldog Applier / Introducer

	Shaft	Working Length	Total Length
9009-395	8 mm	13 3/4" (35 cm)	17 1/3" (44 cm)
9009-393	8 mm	5" (12.7 cm)	9" (23 cm)




SCANLAN® Reliance Bulldog Clamps

Reliable, Surgeon-Controlled Vessel Clamps



2133-03



2135-03

SCANLAN®

Recommended Reliance Bulldog Starter Pack

Recommended	QTY	Cat. No.	Description	Total Length
SCANLAN® V.I.P.™ Sterilization Tray	Choose 1	2133-03	SCANLAN® V.I.P.™ Micro Sterilization Tray With Scan-Guard™ Silicone Mat - Includes Hardware (Holds 3-4 Bulldogs)	1.5" x 3.5" x 1" (3.81 cm x 8.89 cm x 2.54 cm)
		NEW 2135-03	SCANLAN® V.I.P.™ Reliance Bulldog Sterilization Tray With Scan-Guard™ Silicone Mat (Holds 8-10 Bulldogs)	8 1/4" x 3 1/4" x 1 1/2" (21 cm x 8.3 cm x 3.8 cm)
Applier/Introducer	Choose 1	9009-393	Straight Jaws With Ratchet	9" (23 cm)
		9009-395		17 1/3" (44.5 cm)
SCANLAN® Reliance Bulldog Clamps	Up To 4 Of Any Combination	3795-52	40 mm Straight DeBakey Jaws	2 3/4" (7 cm)
		3795-50	25 mm Straight DeBakey Jaws	2" (5 cm)
		NEW 3795-56	15 mm Straight DeBakey Jaws	1 5/8" (4 cm)
		3795-54	7 mm Straight DeBakey Jaws	1 1/4" (3.2 cm)
		3795-51	25 mm Curved DeBakey Jaws	2" (5 cm)
		NEW 3795-57	15 mm Curved DeBakey Jaws	1 5/8" (4 cm)

Firm Sale - May Be Available For Evaluation

For more information or to place an order, please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458
 International: 651-298-0997
 FAX: 651-298-0018
www.scanlaninternational.com



MK-B-RBDC Rev. G