

SCANLAN®


M.D.™ Retractor System

With Radiolucent Scan-Vu™ Arms




Surgical Products Designed and Manufactured
by the Scanlan Family Since 1921

SCANLAN®
M.D.™ Retractor System



Two Rack Sizes
To allow proper patient fit.



Unique Patented Locking System
for a secure, reliable hold of interchangeable arms.

Interchangeable Retractor Arms
Five unique sizes to ensure proper patient fit.



**Scan-Vu™
Radiolucent Arms (Single-Use)**



**M.D.™ Titanium
Retractor Arms (Reusable)**
Designed for light weight, reliable retraction for open surgical procedures.

Durability
Designed with the highest quality materials for reliability & longevity.

SCANLAN® M.D.™ RETRACTOR SYSTEM was developed in cooperation with Francis X. Moga, MD and Ginny Dixon, BSN, RN, CNOR, RNFA

Patent No. 10383613



SCANLAN®
M.D.™ Retractor System



SCANLAN®
M.D.™ Retractor System Size 1 Set 8888-901

Cat. No.	Description	Arm Length	Blade Depth	Blade Width
8008-903	M.D.™ Titanium Rack, 3" (76 mm) Spread			
8008-911*	SCANLAN® Scan-Vu™ Radiolucent Single-use Retractor Arms Size 1 Left and Right (1 pair per sterile package)	1.8" (46 mm)	0.37" (9.5 mm)	0.5" (13 mm)
8008-931R	SCANLAN® M.D.Titanium Retractor Arm Size 1 Right Arm			
8008-931L	SCANLAN® M.D.Titanium Retractor Arm Size 1 Left Arm			



SCANLAN®
M.D.™ Retractor System Size 2 Set 8888-902

Cat. No.	Description	Length	Blade Depth	Blade Width
8008-903	M.D.™ Titanium Rack, 3" (76 mm) Spread			
8008-912*	SCANLAN® Scan-Vu™ Radiolucent Single-use Retractor Arms Size 2 Left and Right (1 pair per sterile package)	2.3" (58 mm)	0.63" (16 mm)	0.9" (23 mm)
8008-932R	SCANLAN® M.D.Titanium Retractor Arm, Size 2 Right Arm			
8008-932L	SCANLAN® M.D.Titanium Retractor Arm, Size 2 Left Arm			



SCANLAN® M.D.™ Retractor System



SCANLAN®
M.D.™ Retractor System Size 3 Set 8888-903

Cat. No.	Description	Length	Blade Depth	Blade Width
8008-905	M.D.™ Titanium Rack, 5" (127 mm) Spread			
8008-913*	SCANLAN® Scan-Vu™ Radiolucent Single-use Retractor Arms Size 3 Left and Right (1 pair per sterile package)	3.0" (76 mm)	0.87" (22 mm)	1.4" (36 mm)
8008-933R	SCANLAN® M.D.Titanium Retractor Arm, Size 3 Right Arm			
8008-933L	SCANLAN® M.D.Titanium Retractor Arm, Size 3 Left Arm			



SCANLAN®
M.D.™ Retractor System Size 4 Set 8888-904

Cat. No.	Description	Length	Blade Depth	Blade Width
8008-905	M.D.™ Titanium Rack, 5" (127 mm) Spread			
8008-914*	SCANLAN® Scan-Vu™ Radiolucent Single-use Retractor Arms, Size 4 Left and Right (1 pair per sterile package)	3.5" (89 mm)	1.14" (29 mm)	1.8" (45 mm)
8008-934R	SCANLAN® M.D.Titanium Retractor Arm, Size 4 Right Arm			
8008-934L	SCANLAN® M.D.Titanium Retractor Arm, Size 4 Left Arm			




SCANLAN®
M.D.™ Retractor System Size 5 Set 8888-905

Cat. No.	Description	Length	Blade Depth	Blade Width
8008-905	M.D.™ Titanium Rack, 5" (127 mm) Spread			
8008-915*	SCANLAN® Scan-Vu™ Radiolucent Single-use Retractor Arms, Size 5 Left and Right (1 pair per sterile package)	4.3" (109 mm)	1.38" (35 mm)	2.6" (67 mm)
8008-935R	SCANLAN® M.D.Titanium Retractor Arm, Size 5 Right Arm			
8008-935L	SCANLAN® M.D.Titanium Retractor Arm, Size 5 Left Arm			



SCANLAN®
M.D.™ Retractor System





SCANLAN®
V.I.P.™ Sterilization Tray

Cat. No.	Description	Dimension
2151-03	Aluminum with Scan-Guard™	17" x 10.5" x 1.88" (43.2 cm x 26.7 cm x 4.8 cm)

New tray designed to accomodate all size variations for quick selection based on patient size.



SCANLAN®
M.D.™ Retractor System



For more information
or to place an order,
please contact your local
Scanlan representative,
or call...

U.S. & Canada: 800-328-9458
International: 651-298-0997
FAX: 651-298-0018
www.scanlaninternational.com

