

SCANLAN®

CAROTID ENDARTERECTOMY INSTRUMENTS



2142-03

SCANLAN ★

Surgical products designed and manufactured by the Scanlan family since 1921
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlangroup.com

INTERNATIONAL

USA

SCANLAN ★

INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921



We are pleased to introduce the SCANLAN® Carotid Endarterectomy Instrument brochure, a collection of new and existing instrument designs that are frequently used in conventional approach carotid endarterectomy procedures.

If you need additional instrument design options, the SCANLAN® catalog or our website (scanlaninternational.com) offers an extensive selection of additional instrument designs.

Even with our extensive design options, we understand the need to accommodate changing technology and unique surgical procedures. To support this, Scanlan International provides customers with in-house design and instrument modification services – where our experienced artisans can customize each instrument to your specifications.

Since 1921, our family has dedicated itself to providing the highest quality surgical instrumentation and fast, friendly service to all of our customers. Our goal is to be the best... and we accomplish this by providing excellence to surgeons, nurses, and ultimately, their patients.


Timothy M. Scanlan
President
Chief Executive Officer




CAROTID ENDARTERECTOMY INSTRUMENTS

[Clamps](#) | [Forceps](#) | [Needle Holders](#) | [Scissors](#) | [Temporary Occlusion Clips](#) | [Miscellaneous](#) | [Trays](#)

SCANLAN® Clamps

Scanlan Endarterectomy Clamps are individually hand-made with the highest grade titanium and stainless steel. They offer a variety of lengths and unique designs to match your specific application ...From the delicate Micro Kitzmiller to the DeBakey Cooley forceps, these clamps offer strong & reliable occlusion.



SCANLAN® Clamps	Cat. No.	Description	Total Length
	2002-235	SCANLAN® Baby Vascular Clamp 45° Angled Cooley Jaws	5 3/4" (14.5 cm)
	4377-08	SCANLAN® DeBakey Cooley Forceps Carmalt Jaws with Cross Hatch Tips	7 1/2" (19 cm)
	5005-60	SCANLAN® Bailey Forceps 30° Angled Jaws Very Fine Angled Serrated Jaws	7" (18 cm)
	5005-62	SCANLAN® Bailey Forceps 60° Angled Jaws Very Fine Angled Serrated Jaws	7" (18 cm)
	5005-64	SCANLAN® Bailey Forceps 90° Angled Jaws Very Fine Angled Serrated Jaws	7" (18 cm)
	3003-614	SCANLAN® Titanium Micro Clamp 60° Angled DeBakey Jaws 23 mm Jaw Length	4 1/8" (10.5 cm)
	3003-628	SCANLAN® Titanium Micro Kitzmiller Clamp Left Angled DeBakey Jaws 23 mm Jaw Length	4 1/4" (11 cm)
	3003-630	SCANLAN® Titanium Micro Kitzmiller Clamp Right Angled DeBakey Jaws 23 mm Jaw Length	4 1/4" (11 cm)



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | [Forceps](#) | Needle Holders | Scissors | Temporary Occlusion Clips | Miscellaneous | Trays

SCANLAN® Forceps

Scanlan® Forceps have been designed to provide secure atraumatic control of tissue. Manufactured of high quality stainless steel and titanium, they provide remarkable balance, precision, control and tactile response on delicate tissue.



SCANLAN® Forceps

	Cat. No.	Description	Total Length
	4004-43	SCANLAN® DeBakey Forceps 1 mm Straight DeBakey Tips	7 3/4" (19.5 cm)
	4004-45	SCANLAN® DeBakey Forceps 2 mm Straight DeBakey Tips	7 3/4" (19.5 cm)
	4004-268	SCANLAN® DeBakey Needle Pulling Tissue Forceps 1 mm Straight Fine Carbon-Inlay with DeBakey Tips Lightweight	7 1/4" (18.5 cm)
	3003-04	SCANLAN® Titanium DeBakey Dissecting Forceps 2 mm Straight DeBakey Tips	7 3/4" (19.5 cm)
	3003-05	SCANLAN® Titanium DeBakey Dissecting Forceps 1 mm Straight DeBakey Tips	7 3/4" (19.5 cm)
	3003-740	SCANLAN® Titanium Gerald DeBakey Forceps 1 mm Straight DeBakey Tips	7" (18 cm)
	3003-924	SCANLAN® LEGACY Micro Forceps 1 mm Straight DIAMOND DUST™ Platform Ring Tips	8 1/4" (21 cm)



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | [Needle Holders](#) | Scissors | Temporary Occlusion Clips | Miscellaneous | Trays

SCANLAN® Needle Holders

SCANLAN® Needle Holders offer a wide selection of styles and designs. Hand - crafted using superior materials designed to become a virtual extension of the surgeon's hand ...providing remarkable balance and precision handling of suture.



SCANLAN® Needle Holders

	Cat. No.	Description	Total Length
	6006-28	SCANLAN® Micro Needle Holder Straight DIAMOND DUST™ Jaws Recommended Suture Size: 5-0 and Smaller	6 1/4" (16 cm)
	6006-82	SCANLAN® Jacobson Micro Needle Holder Straight DIAMOND DUST™ Jaws With Ratchet Recommended Suture Size: 5-0 and Smaller	7 1/4" (18.5 cm)
	6006-120	SCANLAN® Jacobson Micro Needle Holder Straight DIAMOND DUST™ Jaws With Ratchet Recommended Suture Size: 5-0 and Smaller	8 1/4" (21 cm)
	3003-234-1	SCANLAN® Titanium Jacobson Micro Needle Holder Straight Fine DIAMOND DUST™ Jaws With Ratchet Recommended Suture Size: 7-0 and Smaller	8 1/2" (21.5 cm)
	3003-344	SCANLAN® LEGACY Needle Holder Straight DIAMOND DUST™ Jaws With Ratchet Recommended Suture Size: 7-0	8" (20 cm)
	3003-354	SCANLAN® LEGACY Jacobson Micro Needle Holder Straight DIAMOND DUST™ Jaws With Ratchet Recommended Suture Size: 5-0 and Smaller	8" (20 cm)
	3003-434	SCANLAN® Ryder Needle Holder Straight DIAMOND DUST™ Jaws Recommended Suture Size: 5-0 and Smaller	7 1/2" (19 cm)



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | [Scissors](#) | Temporary Occlusion Clips | Miscellaneous | Trays

SCANLAN® SUPER CUT™ Scissors

SUPER CUT™ Scissors incorporate a unique blade design which produces a razor-sharp edge. A wide selection of cutting angles, tips and lengths – facilitates fine, sharp, smooth cuts in vessels that require precision and are difficult to access. Superior balance and tactile feel provide the ultimate in cutting precision and control.



SCANLAN® Scissors

	Cat. No.	Description	Total Length
	7007-211SC	SCANLAN® Super Cut™ Stevens Tenotomy Scissors Curved Blades Tenotomy Tips	7" (18 cm)
	7007-218-3SC	SCANLAN® Super Cut™ Metzenbaum Scissors Curved Blades 2.5 mm Blunt Tips	7" (18 cm)
	7007-252SC	SCANLAN® Super Cut™ Diethrich-Potts Scissors 45° Angled Fine Blades	7" (18 cm)
	7007-320SC	SCANLAN® Super Cut™ DeBakey Scissors 45° Angled Blades Blunt Tips	6 1/4" (16 cm)

CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | Scissors | [Temporary Occlusion Clips](#) | Miscellaneous | Trays

SCANLAN® VASCU-STATT® Single-Use Bulldog Clips

Vascu-Stat® single-use bulldog clamps are a lightweight, cost-effective alternative to metal bulldog clamps. Unlike metal bulldog clamps (clamping pressure weakens after repeated use), Vascu-Stat® single-use bulldog clamps offer a preset tension – ensuring consistent clamping pressure with each use.



SCANLAN® VASCU-STATT® Applier

	Cat. No.	Description	Total Length
	1001-506	SCANLAN® VASCU-STATT® Appliers Spring Style	8 1/2" (21.5 cm)

SCANLAN® VASCU-STATT® Temporary Occlusion Clips

	Cat. No.	Description
	1001-573	SCANLAN® VASCU-STATT® II Plus Midi, Angled With Gentle-Jaw™ Surface Clamping Pressure 75-80 grams 45° (yellow)
	1001-574	SCANLAN® VASCU-STATT® II Plus Maxi, Straight With Gentle-Jaw™ Surface Clamping Pressure 75-80 grams 45° (yellow)



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | Scissors | [Temporary Occlusion Clips](#) | Miscellaneous | Trays

SCANLAN® Reuseable Temporary Occlusion Clips

SCANLAN® Heifetz™ and Yasargil temporary occlusion clips offer surgeons a SMALL alternative to large, bulky metal or disposable clips with similar closing pressure ranges. A wide selection of clip designs provided to match specific surgical applications, such as internal carotid artery.



SCANLAN® Heifetz™ Applier

	Cat. No.	Description	Total Length
	3795-01	SCANLAN® Heifetz™ Appliers Spring Style	7 1/4" (18.5 cm)

SCANLAN® Heifetz™ Temporary Occlusion Clips

	Cat. No.	Description
	3795-18	SCANLAN® Heifetz™ Clip - Firm Pressure Curved 8 mm Jaw Length 1.75 mm Jaw Width
	3795-22	SCANLAN® Heifetz™ Clip - Light Pressure Curved 12 mm Jaw Length 2.25 mm Jaw Width
	3795-28	SCANLAN® Heifetz™ Clip - Firm Pressure Curved 12 mm Jaw Length 2.25 mm Jaw Width



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | Scissors | [Temporary Occlusion Clips](#) | Miscellaneous | Trays

SCANLAN® Heifetz™ & Yasargil Temporary Occlusion Clips

- **Reusable / Durable**
Allowing multiple sterilization cycles, resulting in low cost per procedure
- **Low Profile Design**
Providing excellent visualization
- **Two Color-Coded Ranges of Closing Pressures**
Offering reliable pre-set tensions
- **Specially Designed Appliers**
Providing controlled placement of clips
- **Small Clip Applier Surface**
Reducing suture entanglements



SCANLAN® Yasargil Applying Forceps

	Cat. No.	Description	Total Length
	3786-28	SCANLAN® Yasargil Applying Forceps Bayonet Style	7 1/8" (18.1 cm)

SCANLAN® Yasargil Temporary Occlusion Clips

	Cat. No.	Description
	3785-24	SCANLAN® Yasargil Clip - Light Pressure Slight Curve 9 mm
	3786-24	SCANLAN® Yasargil Clip - Firm Pressure Slight Curve 9 mm
	3786-26	SCANLAN® Yasargil Clip - Firm Pressure Straight 11 mm
	3786-27	SCANLAN® Yasargil Clip - Firm Pressure Slight Curve 11 mm



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | Scissors | Temporary Occlusion Clips | [Miscellaneous](#) | Trays



SCANLAN® Micro Nerve And Vessel Hook

	Cat. No.	Description	Total Length
	3784-32	SCANLAN® Krayenbuhl Micro Nerve and Vessel Hook Small	7 1/4" (18.5 cm)

SCANLAN® Bulldog Shunt Clamp

	Cat. No.	Description	Total Length
	4384-44	SCANLAN® Loftus™ Titanium Bulldog Shunt Clamp 45° Angled 12 Fr.	3 3/4" (9.5 cm)
	4384-46	SCANLAN® Loftus™ Titanium Bulldog Shunt Clamp 60° Angled 12 Fr.	3 3/4" (9.5 cm)
	4385-44	SCANLAN® Loftus™ Titanium Bulldog Shunt Clamp 45° Angled 10 Fr.	3 3/4" (9.5 cm)
	4385-46	SCANLAN® Loftus™ Titanium Bulldog Shunt Clamp 60° Angled 10 Fr.	3 3/4" (9.5 cm)

SCANLAN® Vascular Dilators And Vessel Probes



	Cat. No.	Description	Total Length
	8008-462	SCANLAN® MEMORY Garrett Vascular Dilators & Vessel Probes 3.5 mm	5 1/2" (14 cm)
	8008-464	SCANLAN® MEMORY Garrett Vascular Dilators & Vessel Probes 4 mm	5 1/2" (14 cm)
	8008-466	SCANLAN® MEMORY Garrett Vascular Dilators & Vessel Probes 4.5 mm	5 1/2" (14 cm)
	8008-468	SCANLAN® MEMORY Garrett Vascular Dilators & Vessel Probes 5 mm	5 1/2" (14 cm)



CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | Scissors | Temporary Occlusion Clips | [Miscellaneous](#) | Trays

SCANLAN® Dissectors

	Cat. No.	Description	Total Length
	9009-142	SCANLAN® Ochsner Dissector Angled Delicate Plane Finder Dissector, Angled With Ball	7 1/2" (19 cm)
	9009-144	SCANLAN® Ochsner Dissector Double-Ended Dissectors 1 mm & 2 mm	7 1/8" (18.1 cm)
	9009-146	SCANLAN® Ochsner Dissector Double-Ended Dissectors 3 mm & 4 mm	7 1/8" (18.1 cm)


SCANLAN® Plaque Remover

	Cat. No.	Description	Total Length
	9009-150	SCANLAN® Plaque Remover Curved Fine Tips Ring Handle	6 1/8" 15.5 cm)

SCANLAN® Suction

	Cat. No.	Description	Total Length
	9009-960	SCANLAN® SUNDT™ Graduated Suction System Rotatable Suction Tip Set 3 Fr., 4 Fr., 5 Fr., 6 Fr., 7 Fr., 9 Fr., 11 Fr.	5 1/4" (13.5 cm)

SCANLAN® Vessel Loops

	Cat. No.	Description	Total Length
	1001-76	SCANLAN® SURG-I-LOOP® Maxi Red	16 1/2" (42 cm)
	1001-77	SCANLAN® SURG-I-LOOP® Mini Blue	




CAROTID ENDARTERECTOMY INSTRUMENTS

Clamps | Forceps | Needle Holders | Scissors | Temporary Occlusion Clips | Miscellaneous | [Trays](#)

SCANLAN® V.I.P.™ Metal Sterilization Trays

SCANLAN® V.I.P.™ (Valued Instrument Protection) Sterilization Trays have been designed to offer excellent protection for a variety of surgical instrumentation & devices. Manufactured of strong, surgical grade aluminum, these trays provide reliable sterilization, storage and transportation. Because of their grade durability, they will withstand repeated gas, flash, steam or wrapped sterilization processes.

For added protection, the trays incorporate a SCANLAN® Scan-Guard™ silicone mat and top brackets to securely hold instruments and assure safety.



SCANLAN® V.I.P.™

Metal Sterilization Trays

	Cat. No.	Description	Capacity
	2142-03	Silver-gray color, Scan-Guard™ silicone insert across full width of base, additional tray insert 10" x 15" x 2.6" (25.4 cm x 38 cm x 6.6 cm)	18+ Instruments
	2150-03	Large Capacity Tray Scan-Guard™ silicone brackets & mat 29" x 8 3/8" x 1 1/2" (73.7 cm x 21.3 cm x 3.8 cm)	20 Instruments

For more information or
to place an order,
please contact your local
Scanlan representative,
or call...

U.S. & Canada: 800-328-9458
International: 651-298-0997
FAX: 651-298-0018
www.scanlaninternational.com

