


SCANLAN®
EASY-TAG™
Tracking System



SCANLAN®

EASY-TAG™ Pre-Printed

For even more customer convenience, Scanlan offers pre-printed Easy-Tags.

EASY-TAG™ Pre-Printed

250 tags per spool

2.5" x 4.75" (63.5 mm x 120.65 mm)

Cat. #	Description
1001-921	Blue, Endoscope Cleaned
1001-922	Green, Endo Cleaned Month / Date
1001-923	Red, Rapid Turnover
1001-924	Red, Repair
1001-925	Green, Clean Ready To Use
1001-926	Black, Loaner
1001-927	Red, Biohazard
1001-928	Red, Broken Do Not Use
1001-929	Orange, Tray Is Incomplete
1001-930	Red, Endoscope Repair
1001-931	Green, Process Completed
1001-936	Red, Repair



1001-921

ENDOSCOPE CLEANED BY _____

DATE CLEANED _____

SCOPE SERIAL # _____

EXPIRATION DATE _____

1001-922

JAN • FEB • MAR • APR • MAY • JUN • JUL • AUG • SEP • OCT • NOV • DEC

01 • 02 • 03 • 04 • 05 • 06 • 07

08 • 09 • 10 • 11 • 12 • 13 • 14

15 • 16 • 17 • 18 • 19 • 20 • 21

22 • 23 • 24 • 25 • 26 • 27 • 28

29 • 30 • 31

SCOPE SERIAL # _____ INITIALS _____

1001-923

RAPID
TURNOVER

ROOM # _____ TIME NEEDED _____

SURGEON _____

1001-924

REPAIR NEEDED

Item Broken: _____ Date: _____

Description of Problem: _____

Initial: _____ Critical: Yes or No

1001-925

CLEAN
READY
TO USE

CLEANED BY _____

DATE _____

1001-926

LOANER

VENDOR: _____

DATE NEEDED: _____

OF PANS / TRAYS: _____

FOR DR: _____

FOR ROOM: _____

1001-927

BIOHAZARD

1001-928

BROKEN - DO NOT USE

Date: _____ Reported By: _____

Name: _____ Function: _____

Describe Damage/Problem: _____

1001-929

THIS TRAY IS INCOMPLETE

MISSING: _____

DATE: _____ INITIALS: _____

1001-930

ENDOSCOPE
REPAIR TAG

☐ Blurry ☐ Dark

☐ Cloudy ☐ Bent

☐ Chipped Lens

☐ Light Cord Broken

☐ Wrong Light Cord

☐ Fluid Invasion

☐ Failed Leak Test

☐ Other _____

Completed By: _____

Serial #: _____

Date: _____

1001-931

High Level Disinfection
Process Completed

Ready for Use

Date: _____

Verified by: _____

Serial Number: _____

Exp. Date: _____

1001-936

REPAIR

DATE: _____

ITEM: _____

STAFF NAME: _____

PROBLEM: _____





SCANLAN®

EASY-TAG™ Custom Printing

Customize your EASY-TAG™ Tracking System to meet the unique needs within your hospital. Simply define your message, and we do the rest! Your custom print EASY-TAG™ will be developed and printed to your specific needs.

Your custom message, combined with the unique advantages of the EASY-TAG™ Tracking System, offers a powerful combination! – Providing a cost effective, secure method of tracking instruments, trays, equipment & repairs through cleaning, decontamination and sterilization processes.

Custom EASY-TAG™

2.5" x 4.75" (63.5 mm x 120.65 mm)

Cat. # Description

1001-904C	Black
1001-905C	Red
1001-906C	Yellow
1001-907C	Blue
1001-908C	Green
1001-909C	Pink
1001-910C	Orange

Not made with natural rubber latex
Custom made at time of order. Firm sale.

Custom EASY-TAG™ Small

2.5" x 2.875" (63.5 mm x 73.03 mm)

Cat. # Description


1001-914C	White
1001-915C	Red
1001-916C	Purple
1001-917C	Light Blue
1001-918C	Lime Green
1001-919C	Hot Pink
1001-920C	Orange

Not made with natural rubber latex
Custom made at time of order. Firm sale.



SCANLAN® EASY-TAG™ Pre-Printed Mini

The EASY-TAG™ Mini is Designed to answer the need for a cost-effective, easy, secure method of tracking instruments, trays, equipment & repairs through cleaning, decontamination and sterilization processes.



EASY-TAG™ Pre-Printed Mini

250 Tags Per Box
2.5" x .5" (57.15 mm x 12.7 mm)

Cat. # Description

1001-971	Green, CLEAN
1001-972	Green, LOANER
1001-973	Red, REPAIR
1001-974	Yellow, SHARPEN
1001-975	Orange, BROKEN
1001-976	Blue, RAPID TURNOVER
1001-977	Black, DIRTY



SCANLAN® EASY-TAG™ Plus

Patented flexible/elastic handle technology of the Easy-Tag™ Tracking System, PLUS an adhesive back allowing the user to attach a machine printout, note, or message.

- **Easy to Apply** - the flexible handle loops over endoscope knobs, around instruments, tray handles and carts.
- **Adhesive Back** - allows placement of a machine print out, note or message to the elastic handle
- **Not made with natural rubber latex**
- 500 tags per spool



EASY-TAG™ Plus

Cat. # Description

1001-940	Silver *
----------	----------

* Not intended for sterilization





SCANLAN® EASY-TAG™ Tracking System

A cost-effective, easy, secure method of tracking instruments, trays, equipment & repairs through cleaning, decontamination and sterilization processes.

- **Sterilizable** - EO, Steam and Gas Plasma
- **Versatile Surface** - allows placement of pre-printed labels or can be written on both sides of tag
- **No Adhesives** - quick, easy and secure
- **Easy to Apply** - the flexible handle loops around instruments and tray handles
- **Easy to Remove** - Tear-Away™ handle provides easy separation after use
- **Color Coded** - for fast and easy identification
- **Not made with natural rubber latex**
- 250 Tags per spool

EASY-TAG™



EASY-TAG™ SMALL



EASY-TAG™ Mini



EASY-TAG™

2.5" x 4.75"
(63.5 mm x 120.65 mm)

Cat. #	Description
1001-904	Black
1001-905	Red
1001-906	Yellow
1001-907	Blue
1001-908	Green
1001-909	Pink
1001-910	Orange

EASY-TAG™ Small

2.5" x 2.875"
(63.5 mm x 73.03 mm)

Cat. #	Description
1001-914	White
1001-915	Red
1001-916	Purple
1001-917	Light Blue
1001-918	Lime Green
1001-919	Hot Pink
1001-920	Orange




EASY-TAG™ Mini

250 Tags Per Box
2.5" x .5"
(57.15 mm x 12.7 mm)

Cat. #	Description
1001-964	Black
1001-965	Red
1001-966	Yellow
1001-967	Blue
1001-968	Green
1001-969	Pink
1001-970	Orange



Easy-Tag™ Tracking System Application Guide

Application	Blank Regular	Blank Small	Blank Mini	Pre-Printed	Pre-Printed Mini	Custom Printed	Plus
							
Repair/Broken	✓	✓	✓	1001-924 1001-928 1001-936	1001-973 1001-975 1001-974	✓	
Quick Turnaround	✓	✓	✓	1001-923	1001-976		
Missing Instruments	✓	✓	✓	1001-929		✓	
Loaner Trays	✓	✓	✓	1001-926	1001-972		
Tray/Instrument Identification	✓	✓	✓			✓	
Barcode Placement	✓	✓	✓				
Case Carts	✓	✓	✓	1001-925 1001-927 1001-921 1001-922 1001-925 1001-931		✓	
Endoscope Processing & Hang-Time	✓	✓	✓			✓	✓
Wheel Chairs, IV Poles	✓	✓	✓	✓	1001-971 1001-977	✓	

