


SCANLAN®

LOFTUS™ 2-Level ACDF Distractor



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

SCANLAN® LOFTUS™ 2-Level ACDF Distractor



The Loftus™ 2-Level ACDF Distractor is designed to facilitate a single set-up on two level ACDF procedures. The traditional single level distractor requires the system to be re-positioned and manipulated to accommodate each level. The innovative Loftus™ 2-Level Distractor design overcomes the challenges associated with re-positioning traditional single level distractors.

Additional benefits include:

- **Reduces Operating Room Time** – Incorporates an additional distractor arm permitting a two level set-up to reduce surgical time required for re-positioning and manipulation
- **Cost Effective** – Compatible with existing single level instrumentation and distraction screws
- **Universal** – The shortened bar and center arm enable the distractor to be positioned from either a left or right surgical approach
- **Single Level Manipulation** – Permits independent opening of each cervical space
- **Cleaning and Sterilization** – Distractor arms separate easily for cleaning and sterilization
- **Quality** – Manufactured from high quality stainless steel to provide years of use

Cat. No.	DESCRIPTION
8008-81	Loftus™ 2-Level ACDF Distractor



"I believe this distractor represents a significant advance in patient safety and surgical facility. The patient experiences a safer operation because the old style distractor needs to be removed and replaced to perform a second level surgery. Each removal and manipulation increases the risk of damage to the vital structures in the neck including the esophagus, trachea, carotid artery and jugular vein."

Christopher M. Loftus MD, FACS
Professor and Chairman
Department of Neurosurgery
Temple University School of
Medicine

**For more information
or to place an order,**
please contact your local
Scanlan representative, or call...
U.S. & Canada: 800-328-9458
International: 651-298-0997
FAX: 651-298-0018
www.scanlaninternational.com

