

TRAVELER™

Autoclavable Pouch



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

Traveler™

Autoclavable Pouch

The Scanlan® *Traveler*™ Autoclavable Pouch is designed to provide a fast and reliable method of securing important documents to the outside of containers.

Made of durable material, it is semi-transparent and features a spacious pouch to easily contain large documents. Custom labels may also be applied to the exterior of the pouch to support the identification and tracking of the contents.

Features Include...

- Sterilizable – Autoclave sterilization tested
- Sized to Fit – Holds instrument count sheet with only minimal folding
- Contents Secure – Fastener and dual hanging hole ensure contents of pouch do not fall out
- Easy On – Adjustable fastener makes it easy to put on and take off, no scissors and no adhesives to leave residue behind
- Durable – Strong 3 mil material
- Large Surface Area – Perfect for placing a bar code label or to write on
- Versatile – Autoclavable pouch is great for housing various documents required to go through the sterilization process with the container

Order Information

| Catalog Number | Description |
|----------------|---|
| 1001-950 | Traveler™ Autoclavable Pouch 100 pouches (6.5" x 6.5") and 100 fasteners/box |
| 1001-951 | Traveler™ Autoclavable Pouch 100 pouches (6.5" x 6.5") (fasteners not included) |

* Patent pending



Insert count sheet



Labels may be attached (not included)



Insert the fastener through the autoclavable pouch



Place around the container handle and adjust the fastener



Autoclavable pouch and container

For more information or to place an order,
please contact your local Scanlan representative,

or call... 800-328-9458
FAX: 651-298-0018
www.scanlaninternational.com



INTERNATIONAL