

SURG-I-BAND®

Color Coding Data Matrix Barcode



SURG-I-BAND® Color Coding Data Matrix Barcode

Offering Precise Instrument / Device Tracking

Surg-I-Band® Data Matrix Color Coding Barcode offers hospitals an excellent system to organize and track personal protective gear, surgical instrumentation and devices. Combining the benefits of Surg-I-Band® color coding with the technology of bar coding, Surg-I-Band® Data Matrix provides a fast, easy and precise organizational system.

- **Economy** – Each spool contains 200 unique data matrix barcodes
- **14 colors** – For fast and easy identification
- **Safety** – Clinically tested for steam, ETO and Immediate Use Steam Sterilization (IUSS)
- **Trouble Free** – Easy application and removal without damaging instrument
- **Versatility** – Can be used on most sizes and designs of instruments
- Compatible with most instrument tracking software and high-density data matrix ECC 200 scanners
- May be used in combination with Surg-I-Band® Dispenser, Remover and Carousel
- 1/4" Width
- Not made with natural rubber latex




























SCANLAN®

SURG-I-BAND® Data Matrix Barcode Pre-Cut Assorted Sizes

Designed for tracking and barcoding larger diameter instruments.

Pre-cut label sheet (8 1/2" x 11") - 1/4" width, 132 barcodes per sheet

- 72 each of 1.75"
- 27 each of 3"
- 26 each of 4"
- 7 each of 5"

Description	Censitrac™ Compatible	Assorted Sizes Sheets
Red	 	1001-1098 1001-1112
Yellow	 	1001-1099 1001-1113
Orange	 	1001-1100 1001-1114
Blue	 	1001-1101 1001-1115
Green	 	1001-1102 1001-1116
Brown	 	1001-1103 1001-1117
White	 	1001-1104 1001-1118
Black	 	1001-1105 1001-1119
Glo-Orange	 	1001-1106 1001-1120
Purple	 	1001-1107 1001-1121
Lime Green	 	1001-1108 1001-1122
Grey	 	1001-1109 1001-1123
Pink	 	1001-1110 1001-1124
Light Blue	 	1001-1111 1001-1125

