

# SINGLE-USE & INSTRUMENT CARE

Products



SCANLAN ★

INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*



SCANLAN®  
**SURG-I-BAND®**  
Color Coding



Dear Customer,


Since 1921, Scanlan International has supplied hospital operating rooms with the highest quality surgical instrumentation and products. Since inventing Surg-I-Band® color coding (our first instrument care product), Scanlan has worked hand-in-hand with surgeons and perioperative nurses around the world to develop new and innovative products.

Today, Scanlan offers more than 25 single-use and instrument care products – each reflecting the Scanlan commitment to quality through superior designs, the highest quality medical-grade materials, rigorous evaluation and testing.

Our goal is simple – provide the finest surgical products to assist the world’s surgeons, nurses and ultimately... their patients.

Sincerely,

Timothy M. Scanlan  
President  
Chief Executive Officer



Original Surg-I-Band® (1/4") and narrow Surg-I-Band® (1/8") color coding, available in 14 bright colors, offers a multitude of color coding combinations to identify instruments by surgeon, set, specialty and suite. In addition, Surg-I-Band® color coding is long lasting – providing greater value by saving time and labor. Not made with natural rubber latex. Tested for steam, ETO and Immediate Use Steam Sterilization (IUSS).

| 1/4"<br>Original<br>Cat. No. | Description |  | 1/8"<br>Narrow<br>Cat. No. | Description |  |
|------------------------------|-------------|--|----------------------------|-------------|--|
| 1001-64                      | Red         |  | 1001-630                   | Red         |  |
| 1001-65                      | Yellow      |  | 1001-631                   | Yellow      |  |
| 1001-66                      | Orange      |  | 1001-632                   | Orange      |  |
| 1001-67                      | Blue        |  | 1001-633                   | Blue        |  |
| 1001-68                      | Green       |  | 1001-634                   | Green       |  |
| 1001-69                      | Brown       |  | 1001-635                   | Brown       |  |
| 1001-70                      | White       |  | 1001-636                   | White       |  |
| 1001-71                      | Black       |  | 1001-637                   | Black       |  |
| 1001-94                      | Glo-Orange  |  | 1001-638                   | Glo-Orange  |  |
| 1001-95                      | Purple      |  | 1001-639                   | Purple      |  |
| 1001-96                      | Lime Green  |  | 1001-640                   | Lime Green  |  |
| 1001-97                      | Grey        |  | 1001-641                   | Grey        |  |
| 1001-98                      | Pink        |  | 1001-642                   | Pink        |  |
| 1001-99                      | Light Blue  |  | 1001-643                   | Light Blue  |  |





## SCANLAN® SURG-I-BAND® Designs On Spools

- **40 designs** – numbers, symbols and abbreviations
- **Multiple colors**
- **Durable** – resistant to flaking, chipping, peeling and cracking through multiple sterilization cycles
- **Sterilizable** – Tested for steam, ETO and Immediate Use Steam Sterilization (IUSS)
- **1/4" Width**
- **Not made with natural rubber latex**



### Red Designs

| Cat. No.  | Description |             |
|-----------|-------------|-------------|
| 1001-1000 | Zero        | 000000000   |
| 1001-1001 | One         | 11111111111 |
| 1001-1002 | Two         | 222222222   |
| 1001-1003 | Three       | 333333333   |
| 1001-1004 | Four        | 444444444   |
| 1001-1005 | Five        | 555555555   |
| 1001-1006 | Six         | 666666666   |
| 1001-1007 | Seven       | 777777777   |
| 1001-1008 | Eight       | 888888888   |
| 1001-1009 | Nine        | 999999999   |
| 1001-1010 | Hearts      | ♥♥♥♥♥♥♥♥♥   |
| 1001-1011 | Bones       | ☠☠☠☠☠☠☠☠☠   |
| 1001-1012 | Crosses     | ✚✚✚✚✚✚✚✚✚   |
| 1001-1013 | Eyes        | 👁👁👁👁👁👁👁👁    |

### Pink Designs

| Cat. No.  | Description  |                         |
|-----------|--------------|-------------------------|
| 1001-1042 | OB           | OB OB OB OB OB OB OB    |
| 1001-1043 | EYE          | EYE EYE EYE EYE EYE EYE |
| 1001-1044 | GEN          | GEN GEN GEN GEN GEN GEN |
| 1001-1045 | GYN          | GYN GYN GYN GYN GYN GYN |
| 1001-1046 | BAR          | BAR BAR BAR BAR BAR BAR |
| 1001-1047 | URO          | URO URO URO URO URO URO |
| 1001-1048 | NEU          | NEU NEU NEU NEU NEU NEU |
| 1001-1049 | L&D          | L&D L&D L&D L&D L&D L&D |
| 1001-1050 | Feet         | 👣👣👣👣👣👣👣👣                |
| 1001-1051 | Baby Bottle  | 🍼🍼🍼🍼🍼🍼🍼🍼                |
| 1001-1052 | Bones        | ☠☠☠☠☠☠☠☠☠               |
| 1001-1053 | Smiley Faces | 😊😊😊😊😊😊😊😊😊               |
| 1001-1054 | Hearts       | ♥♥♥♥♥♥♥♥♥               |
| 1001-1055 | Hands        | 👐👐👐👐👐👐👐👐👐               |

### Black Designs

| Cat. No.  | Description |             |
|-----------|-------------|-------------|
| 1001-1014 | Zero        | 000000000   |
| 1001-1015 | One         | 11111111111 |
| 1001-1016 | Two         | 222222222   |
| 1001-1017 | Three       | 333333333   |
| 1001-1018 | Four        | 444444444   |
| 1001-1019 | Five        | 555555555   |
| 1001-1020 | Six         | 666666666   |
| 1001-1021 | Seven       | 777777777   |
| 1001-1022 | Eight       | 888888888   |
| 1001-1023 | Nine        | 999999999   |
| 1001-1024 | Hearts      | ♥♥♥♥♥♥♥♥♥   |
| 1001-1025 | Bones       | ☠☠☠☠☠☠☠☠☠   |
| 1001-1026 | Crosses     | ✚✚✚✚✚✚✚✚✚   |
| 1001-1027 | Eyes        | 👁👁👁👁👁👁👁👁    |

### Green Designs

| Cat. No.  | Description |                               |
|-----------|-------------|-------------------------------|
| 1001-1056 | ETO         | ETO ETO ETO ETO ETO ETO       |
| 1001-1057 | ENDO        | ENDO ENDO ENDO ENDO ENDO ENDO |
| 1001-1058 | GI          | GI GI GI GI GI GI GI GI       |
| 1001-1059 | CV          | CV CV CV CV CV CV CV CV       |
| 1001-1060 | ER          | ER ER ER ER ER ER ER ER       |
| 1001-1061 | ASC         | ASC ASC ASC ASC ASC ASC ASC   |
| 1001-1062 | CCU         | CCU CCU CCU CCU CCU CCU CCU   |
| 1001-1063 | DEN         | DEN DEN DEN DEN DEN DEN DEN   |
| 1001-1064 | ENT         | ENT ENT ENT ENT ENT ENT ENT   |
| 1001-1065 | ICU         | ICU ICU ICU ICU ICU ICU ICU   |
| 1001-1066 | SPD         | SPD SPD SPD SPD SPD SPD SPD   |
| 1001-1067 | Triangles   | ▲▲▲▲▲▲▲▲▲                     |
| 1001-1068 | Circles     | ●●●●●●●●●                     |
| 1001-1069 | Stars       | ★★ ★★ ★★ ★★ ★★ ★★             |

### Blue Designs

| Cat. No.  | Description  |                         |
|-----------|--------------|-------------------------|
| 1001-1028 | OB           | OB OB OB OB OB OB OB    |
| 1001-1029 | EYE          | EYE EYE EYE EYE EYE EYE |
| 1001-1030 | GEN          | GEN GEN GEN GEN GEN GEN |
| 1001-1031 | GYN          | GYN GYN GYN GYN GYN GYN |
| 1001-1032 | BAR          | BAR BAR BAR BAR BAR BAR |
| 1001-1033 | URO          | URO URO URO URO URO URO |
| 1001-1034 | NEU          | NEU NEU NEU NEU NEU NEU |
| 1001-1035 | L&D          | L&D L&D L&D L&D L&D L&D |
| 1001-1036 | Feet         | 👣👣👣👣👣👣👣👣                |
| 1001-1037 | Baby Bottle  | 🍼🍼🍼🍼🍼🍼🍼🍼                |
| 1001-1038 | Bones        | ☠☠☠☠☠☠☠☠☠               |
| 1001-1039 | Smiley Faces | 😊😊😊😊😊😊😊😊😊               |
| 1001-1040 | Hearts       | ♥♥♥♥♥♥♥♥♥               |
| 1001-1041 | Hands        | 👐👐👐👐👐👐👐👐👐               |

### Purple Designs

| Cat. No.  | Description |                               |
|-----------|-------------|-------------------------------|
| 1001-1070 | ETO         | ETO ETO ETO ETO ETO ETO       |
| 1001-1071 | ENDO        | ENDO ENDO ENDO ENDO ENDO ENDO |
| 1001-1072 | GI          | GI GI GI GI GI GI GI GI       |
| 1001-1073 | CV          | CV CV CV CV CV CV CV CV       |
| 1001-1074 | ER          | ER ER ER ER ER ER ER ER       |
| 1001-1075 | ASC         | ASC ASC ASC ASC ASC ASC ASC   |
| 1001-1076 | CCU         | CCU CCU CCU CCU CCU CCU CCU   |
| 1001-1077 | DEN         | DEN DEN DEN DEN DEN DEN DEN   |
| 1001-1078 | ENT         | ENT ENT ENT ENT ENT ENT ENT   |
| 1001-1079 | ICU         | ICU ICU ICU ICU ICU ICU ICU   |
| 1001-1080 | SPD         | SPD SPD SPD SPD SPD SPD SPD   |
| 1001-1081 | Triangles   | ▲▲▲▲▲▲▲▲▲                     |
| 1001-1082 | Circles     | ●●●●●●●●●                     |
| 1001-1083 | Stars       | ★ ★★ ★★ ★★ ★★ ★★              |

## SCANLAN® SURG-I-BAND® Color Coding Labels

- **Three patterns** – polka-dots, checkered and Jammie-Stripes™
- **14 solid colors**
- **Durable** – resistant to flaking, chipping, peeling and cracking through multiple sterilization cycles
- **Sterilizable** – Tested for steam, ETO and Immediate Use Steam Sterilization (IUSS)
- **Pre-cut label sheet** (8 1/2" x 11") – 1/4" width, 374" per sheet
  - 230 ea. of 1.0"
  - 40 ea. of 1.5"
  - 12 ea. of 2.0"
  - 12 ea. of 5.0"
- **Not made with natural rubber latex**

### Solid Labels

| Cat. No. | Description |  |
|----------|-------------|--|
| 1001-800 | Red         |  |
| 1001-801 | Yellow      |  |
| 1001-802 | Orange      |  |
| 1001-803 | Blue        |  |
| 1001-804 | Green       |  |
| 1001-805 | Brown       |  |
| 1001-806 | White       |  |
| 1001-807 | Black       |  |
| 1001-808 | Glo-Orange  |  |
| 1001-809 | Purple      |  |
| 1001-810 | Lime Green  |  |
| 1001-811 | Grey        |  |
| 1001-812 | Pink        |  |
| 1001-813 | Light Blue  |  |

### Polka-Dot Labels

| Cat. No. | Description  |  |
|----------|--------------|--|
| 1001-814 | Red/White    |  |
| 1001-815 | Yellow/White |  |
| 1001-816 | Green/White  |  |
| 1001-817 | Black/White  |  |

### Checked Labels

| Cat. No. | Description  |  |
|----------|--------------|--|
| 1001-818 | Red/White    |  |
| 1001-819 | Yellow/White |  |
| 1001-820 | Green/White  |  |
| 1001-821 | Black/White  |  |

### Jammie-Stripes Labels

| Cat. No. | Description  |  |
|----------|--------------|--|
| 1001-822 | Red/White    |  |
| 1001-823 | Yellow/White |  |
| 1001-824 | Green/White  |  |
| 1001-825 | Black/White  |  |



SCANLAN®  
**SURG-I-BAND®**  
Color Coding Data Matrix Barcode  
Regular and Censitrac™ Compatible

**Offering Precise Instrument / Device Tracking**

Surg-I-Band® Data Matrix Color Coding Barcode offers hospitals an excellent system to organize and track surgical instrumentation and devices. Combining the benefits of Surg-I-Band® color coding with the technology of bar coding, Surg-I-Band® Data Matrix provides a fast, easy and precise organizational system.

**Advantages Include:**

- **14 colors** – For fast and easy identification
- **Safety** – Tested for steam, ETO and Immediate Use Steam Sterilization (IUSS)
- **Trouble Free** – Easy application and removal without damaging instrument
- **Versatility** – Can be used on most sizes and designs of instruments
- Compatible with instrument tracking software and high-density data matrix ECC 200 scanners
- May be used in combination with Surg-I-Band® Dispenser, Remover and Carousel
- 1/4" Width
- Not made with natural rubber latex

**Available in:**

- **Spools** – contains 200 unique data matrix barcodes
- **Pre-Cut Label Sheets** – (8 1/2" x 11") - 1/4" width, 132 unique data matrix barcodes per sheet
  - 72 each of 1.75"
  - 27 each of 3"
  - 26 each of 4"
  - 7 each of 5"



| Description |                                                                                     | Data Matrix                                                                         | Censitrac™ Compatible | Assorted Sizes Labels<br>(1 Sheet, 1/4" Width) |           |
|-------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------|------------------------------------------------|-----------|
| Red         |    |    | 1001-1084             | 1001-1098                                      | 1001-1112 |
| Yellow      |    |    | 1001-1085             | 1001-1099                                      | 1001-1113 |
| Orange      |    |    | 1001-1086             | 1001-1100                                      | 1001-1114 |
| Blue        |   |   | 1001-1087             | 1001-1101                                      | 1001-1115 |
| Green       |  |  | 1001-1088             | 1001-1102                                      | 1001-1116 |
| Brown       |  |  | 1001-1089             | 1001-1103                                      | 1001-1117 |
| White       |  |                                                                                     | 1001-1090             | 1001-1104                                      | 1001-1118 |
| Black       |  |  | 1001-1091             | 1001-1105                                      | 1001-1119 |
| Glo-Orange  |  |  | 1001-1092             | 1001-1106                                      | 1001-1120 |
| Purple      |  |  | 1001-1093             | 1001-1107                                      | 1001-1121 |
| Lime Green  |  |  | 1001-1094             | 1001-1108                                      | 1001-1122 |
| Grey        |  |  | 1001-1095             | 1001-1109                                      | 1001-1123 |
| Pink        |  |  | 1001-1096             | 1001-1110                                      | 1001-1124 |
| Light Blue  |  |  | 1001-1097             | 1001-1111                                      | 1001-1125 |



SCANLAN®  
**SURG-I-BAND®**  
Accessories

SCANLAN®

**SURG-I-BAND® Carousel**

- Efficient – provides fast, easy application of Surg-I-Band® color coding
- Excellent organization – holds 12 Surg-I-Band® dispensers
- Rotating top – provides easy access to each color
- Storage – conveniently holds Scanlan® Tip Guard™ instrument protectors
- Easy to use – Surg-I-Band® spools snap on securely and can be removed quickly



| Cat. No. | Description                                                                         |
|----------|-------------------------------------------------------------------------------------|
| 1001-59  | Carousel Only                                                                       |
| 1001-61  | Carousel with 12 Dispensers<br>(Does not include Surg-I-Band® spools or Tip-Guard™) |

SCANLAN®

**SURG-I-BAND® Dispenser**



- Easy to use – Surg-I-Band® color coding spools snap on securely and can be removed quickly
- Precise cut – Dispenser eliminates the need for awkward scissors or blades when using Surg-I-Band® color coding
- Protective cover - guards against damage to blade and color coding
- Compact – provides a small, lightweight design
- Reusable – ensuring greater economy

| Cat. No. | Description                                                      |
|----------|------------------------------------------------------------------|
| 1001-60  | Surg-I-Band® Dispenser<br>(Does not include Surg-I-Band® spools) |

SCANLAN®

**SURG-I-BAND® Remover**



- Saves time – allows quick removal of Surg-I-Band® color coding
- Safe – eliminates use of razors or other sharp instruments that may damage the surface of instruments or cause injury to personnel
- Reusable – removes hundreds of bands

| Cat. No. | Description          |
|----------|----------------------|
| 1001-49  | Surg-I-Band® Remover |





## SCANLAN® EASY-TAG™ Tracking & Tagging System



A cost-effective, easy, secure method of tracking instruments, trays, equipment & repairs through cleaning, decontamination and sterilization processes.

- **Sterilizable** - ETO, Steam and Gas Plasma
- **Versatile Surface** - allows placement of pre-printed labels or can be written on both sides of tag
- **No Adhesives** - quick, easy and secure
- **Easy to Apply** - the flexible handle loops around instruments and tray handles
- **Easy to Remove** - Tear-Away™ handle provides easy separation after use
- **Color Coded** - for fast and easy identification
- Not made with natural rubber latex
- Approximately 250 Tags per spool



## SCANLAN® EASY-TAG™ Tracking & Tagging System

### SCANLAN® EASY-TAG™ Tracking System

| Cat. No. | Description       |
|----------|-------------------|
| 1001-904 | EASY-TAG™, Black  |
| 1001-905 | EASY-TAG™, Red    |
| 1001-906 | EASY-TAG™, Yellow |
| 1001-907 | EASY-TAG™, Blue   |
| 1001-908 | EASY-TAG™, Green  |
| 1001-909 | EASY-TAG™, Pink   |
| 1001-910 | EASY-TAG™, Orange |



2.5" x 4.75"  
(63.5 mm x 120.65 mm)

### SCANLAN® EASY-TAG™ Small Tracking System


| Cat. No. | Description           |
|----------|-----------------------|
| 1001-914 | EASY-TAG™, White      |
| 1001-915 | EASY-TAG™, Red        |
| 1001-916 | EASY-TAG™, Purple     |
| 1001-917 | EASY-TAG™, Light Blue |
| 1001-918 | EASY-TAG™, Lime Green |
| 1001-919 | EASY-TAG™, Hot Pink   |
| 1001-920 | EASY-TAG™, Orange     |



2.5" x 2.875"  
(63.5 mm x 73.03 mm)

### SCANLAN® EASY-TAG™ Mini Tracking System

| Cat. No. | Description       |
|----------|-------------------|
| 1001-964 | EASY-TAG™, Black  |
| 1001-965 | EASY-TAG™, Red    |
| 1001-966 | EASY-TAG™, Yellow |
| 1001-967 | EASY-TAG™, Blue   |
| 1001-968 | EASY-TAG™, Green  |
| 1001-969 | EASY-TAG™, Pink   |
| 1001-970 | EASY-TAG™, Orange |



0.5" x 2.25"  
(12.7 mm x 57.15 mm)

Approximately 250 Tags per box



# SCANLAN® EASY-TAG™ Tracking & Tagging System

## SCANLAN® EASY-TAG™ Custom Printing – Available in Multiple Languages

Customize your EASY-TAG™ Tracking System to meet the unique needs within your hospital. Simply define your message, and we do the rest! Your custom print EASY-TAG™ will be developed and printed to your specific needs.

Your custom message, combined with the unique advantages of the EASY-TAG™ Tracking System, offers a powerful combination – providing a cost-effective, secure method of tracking instruments, trays, equipment & repairs through cleaning, decontamination and sterilization processes.



| Cat. No.  | Description       |
|-----------|-------------------|
| 1001-904C | EASY-TAG™, Black  |
| 1001-905C | EASY-TAG™, Red    |
| 1001-906C | EASY-TAG™, Yellow |
| 1001-907C | EASY-TAG™, Blue   |
| 1001-908C | EASY-TAG™, Green  |
| 1001-909C | EASY-TAG™, Pink   |
| 1001-910C | EASY-TAG™, Orange |

Not made with natural rubber latex  
Spool of 250 Tags

Custom made at time of order. Firm sale.

## SCANLAN® EASY-TAG™ Small with Custom Printing

| Cat. No.  | Description           |
|-----------|-----------------------|
| 1001-914C | EASY-TAG™, White      |
| 1001-915C | EASY-TAG™, Red        |
| 1001-916C | EASY-TAG™, Purple     |
| 1001-917C | EASY-TAG™, Light Blue |
| 1001-918C | EASY-TAG™, Lime Green |
| 1001-919C | EASY-TAG™, Hot Pink   |
| 1001-920C | EASY-TAG™, Orange     |

Not made with natural rubber latex  
Spool of 250 Tags

Custom made at time of order. Firm sale.



# SCANLAN® EASY-TAG™ Tracking & Tagging System

## SCANLAN® EASY-TAG™ Pre-Printed

| Cat. No. | Description                                           |
|----------|-------------------------------------------------------|
| 1001-921 | EASY-TAG™, Blue<br>Endoscope Cleaned                  |
| 1001-922 | EASY-TAG™, Green<br>Endoscope Cleaned<br>Month / Date |
| 1001-923 | EASY-TAG™, Red<br>Rapid Turnover                      |
| 1001-924 | EASY-TAG™, Red<br>Repair                              |
| 1001-925 | EASY-TAG™, Green<br>Clean Ready To Use                |
| 1001-926 | EASY-TAG™, Black<br>Loaner                            |
| 1001-927 | EASY-TAG™, Red<br>Biohazard                           |

Not made with natural rubber latex  
Spool of 250 Tags



01

ENDOSCOPE CLEANED BY \_\_\_\_\_  
DATE CLEANED \_\_\_\_\_  
SCOPE SERIAL # \_\_\_\_\_  
EXPIRATION DATE \_\_\_\_\_

02

JAN • FEB • MAR • APR • MAY • JUN • JUL • AUG • SEP • OCT • NOV • DEC  
 01 • 02 • 03 • 04 • 05 • 06 • 07  
 08 • 09 • 10 • 11 • 12 • 13 • 14  
 15 • 16 • 17 • 18 • 19 • 20 • 21  
 22 • 23 • 24 • 25 • 26 • 27 • 28  
 29 • 30 • 31  
 SCOPE SERIAL # \_\_\_\_\_ INITIALS \_\_\_\_\_

03

**RAPID  
TURNOVER**  
 ROOM # \_\_\_\_\_ TIME NEEDED \_\_\_\_\_  
 SURGEON \_\_\_\_\_

04

Item Broken: \_\_\_\_\_ Date: \_\_\_\_\_  
 Description of Problem: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Initials: \_\_\_\_\_ Critical: Yes or No

05

**CLEAN  
READY  
TO USE**  
 CLEANED BY \_\_\_\_\_  
 DATE \_\_\_\_\_

06

**LOANER**  
 VENDOR: \_\_\_\_\_  
 DATE NEEDED: \_\_\_\_\_  
 # OF PANS / TRAYS: \_\_\_\_\_  
 FOR DR: \_\_\_\_\_  
 FOR ROOM: \_\_\_\_\_

07

**BIOHAZARD**





## SCANLAN® EASY-TAG™ Plus Tracking & Tagging System

Patented flexible/elastic handle technology of the Easy-Tag Tracking System, plus an adhesive back allowing the user to attach a machine printout, note, or message.

- **Easy to Apply** - the flexible handle loops over endoscope knobs, around instruments, tray handles and carts.
- **Adhesive Back** - allows placement of a machine print out, note or message to the elastic handle
- **Not made with natural rubber latex**
- **Approximately 500 tags per spool**



| Cat. No.               | Description            |
|------------------------|------------------------|
| <b>NEW</b><br>1001-940 | Easy-Tag Plus, Silver* |

\* Not intended for sterilization



## SCANLAN® TIP-GUARD™ Instrument Protectors

Ensuring safe, cost-effective protection of blade tips and edges for all surgical instrumentation.

- **Economical** - provides long-term use and is cost-effective
- **Superior protection** - safeguards delicate tips and edges of surgical instrumentation
- **Tested safe** - clinically tested for steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization (sterility report available upon request)
- **Versatile** - conforms to various shapes and angles
- **Design variety** - available in colored, tinted, clear and vented
- **Radiopaque**
- **Not made with natural rubber latex**



|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
|---------|-------------------|-------------------|--------------------|--------------------|--------------------|-----------------|-------------------|------------------|-------------------|--------------------------------------------------------------------------------------|
| Colored |                   |                   |                    |                    |                    |                 |                   |                  |                   | <b>1001-84 Variety Pack</b><br>9 assorted color-coded sizes, 100 pieces per package  |
|         | 1001-89<br>Black  | 1001-88<br>Brown  | 1001-87<br>Orange  | 1001-91<br>Maroon  | 1001-90<br>Yellow  | 1001-85<br>Red  | 1001-86<br>Green  | 1001-93<br>Blue  | 1001-92<br>White  |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
| Tinted  |                   |                   |                    |                    |                    |                 |                   |                  |                   | <b>1001-300 Variety Pack</b><br>9 assorted color-coded sizes, 100 pieces per package |
|         | 1001-305<br>Black | 1001-304<br>Brown | 1001-303<br>Orange | 1001-307<br>Maroon | 1001-306<br>Yellow | 1001-301<br>Red | 1001-302<br>Green | 1001-309<br>Blue | 1001-308<br>White |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
| Clear   |                   |                   |                    |                    |                    |                 |                   |                  |                   | <b>1001-200 Variety Pack</b><br>9 assorted sizes, 100 pieces per package             |
|         | 1001-205          | 1001-204          | 1001-203           | 1001-207           | 1001-206           | 1001-201        | 1001-202          | 1001-209         | 1001-208          |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
| Vented  |                   |                   |                    |                    |                    |                 |                   |                  |                   | <b>1001-400 Variety Pack</b><br>9 assorted color-coded sizes, 100 pieces per package |
|         | 1001-405<br>Black | 1001-404<br>Brown | 1001-403<br>Orange | 1001-407<br>Maroon | 1001-406<br>Yellow | 1001-401<br>Red | 1001-402<br>Green | 1001-409<br>Blue | 1001-408<br>White |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |
|         |                   |                   |                    |                    |                    |                 |                   |                  |                   |                                                                                      |



SCANLAN®

V.I.P™ Sterilization Trays (Valued Instrument Protection)

Metal

SCANLAN® instrument sterilization trays assist in the sterilization, storage, and transport of surgical instruments. All trays will withstand repeated ETO, Immediate Use Steam Sterilization (IUSS), steam or wrapped sterilization and include silicone mats to securely hold instruments and assure safety.

SCANLAN®


V.I.P™ Micro Sterilization Trays

Aluminum, Scan-Guard™ Silicone Inserts

| Cat. No.              | Description                                                         | Dimensions                                            | Capacity |
|-----------------------|---------------------------------------------------------------------|-------------------------------------------------------|----------|
| 2133-03               | Micro sterilization tray with Scan-Guard™ silicone mat and hardware | 1.5" x 3.5" x 1"<br>(3.81 cm x 8.89 cm x 2.54 cm)     | 3-4      |
| <b>NEW</b><br>2135-03 | Scan-Guard™ Silicone Mat                                            | 8 1/4" x 3 1/4" x 1 1/2"<br>(21 cm x 8.3 cm x 3.8 cm) | 8-10     |



(Surgical parts not included)  
Firm Sale – May be available for evaluation



SCANLAN®

V.I.P™ Metal Sterilization Trays

Aluminum, Scan-Guard™ Silicone Inserts, Detachable Lids

| Cat. No.              | Description                                                                                                                                       | Dimensions                                            | Capacity                                                                                     |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------|----------------------------------------------------------------------------------------------|
| 2145-03               | VATS / MIS Tray with Scan-Guard™ brackets<br><i>Firm Sale – May be available for evaluation</i>                                                   | 10.5" x 21" x 2.5"<br>(26 cm x 53 cm x 6.4 cm)        | Up to 12 VATS / MIS Instruments                                                              |
| 2146-03               | Scan-Guard™ silicone bracket with additional tray insert with full Scan-Guard™ silicone mat<br><i>Firm Sale – May be available for evaluation</i> | 10.25" x 21.75" x 4"<br>(26 cm x 55.2 cm x 10 cm)     | 20+ VATS / MIS Instruments                                                                   |
| <b>NEW</b><br>2150-03 | Large capacity tray; Scan-Guard™ silicone brackets & mat<br><i>Firm Sale – May be available for evaluation</i>                                    | 29" x 8 3/8" x 1 1/2"<br>(73.7 cm x 21.3 cm x 3.8 cm) | Holds the entire SCANLAN® Tunneling System<br>• 8 VATS / MIS Instruments<br>• 20 Instruments |

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible

*Firm Sale – May be available for evaluation*



SCANLAN®

V.I.P™ Sterilization Trays (Valued Instrument Protection)

Metal



SCANLAN®

V.I.P™ Metal Sterilization Trays

Aluminum, Scan-Guard™ Silicone Inserts, Detachable Lids

| Cat. No. | Description                                                                                         | Dimensions                                       | Capacity                           |
|----------|-----------------------------------------------------------------------------------------------------|--------------------------------------------------|------------------------------------|
| 2138-03  | Scan-Guard™ Silicone brackets and mat inserts                                                       | 10" x 10" x 1.5"<br>(25.4 cm x 25.4 cm x 3.8 cm) | 8-10                               |
| 2139-03  | Full Scan-Guard™ silicone mat on base, mat insert and brackets on cover to secure small instruments | 10" x 15" x 1.5"<br>(25.4 cm x 38 cm x 3.8 cm)   | 20+ or<br>7 VATS / MIS Instruments |
| 2140-03  | Scan-Guard™ Silicone brackets and mat inserts                                                       | 10" x 15" x 1.5"<br>(25.4 cm x 38 cm x 3.8 cm)   | 10-12                              |
| 2141-03  | Base with accessory compartment, additional tray insert, with Scan-Guard™ Silicone Inserts          | 10" x 15" x 2.6"<br>(25.4 cm x 38 cm x 6.6 cm)   | 15+                                |
| 2142-03  | Scan-Guard™ silicone insert across full width of base, additional tray insert                       | 10" x 15" x 2.6"<br>(25.4 cm x 38 cm x 6.6 cm)   | 18+                                |


All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible

*Firm Sale – May be available for evaluation*





SCANLAN®  
**V.I.P™ Sterilization Trays** (Valued Instrument Protection)  
Plastic



SCANLAN®  
**V.I.P™ Plastic Sterilization Trays**  
Scan-Guard™ Silicone Inserts, Detachable Lids


| Cat. No. | Description                                                                                        | Capacity |
|----------|----------------------------------------------------------------------------------------------------|----------|
| 2083-01  | 6" x 2" x 0.7" (15 cm x 5.1 cm x 1.8 cm)                                                           | 1-2      |
| 2083-02  | 6" x 3" x 0.7" (15 cm x 7.6 cm x 1.8 cm)                                                           | 3        |
| 2083-03  | 7.8" x 3" x 0.7" (19.8 cm x 7.6 cm x 1.8 cm)                                                       | 3        |
| 2083-05  | 9.6" x 5.5" x 0.7" (24.4 cm x 14 cm x 1.8 cm)                                                      | 5-10     |
| 2081-11  | 10" x 6" x 0.7" (25.4 cm x 15 cm x 1.8 cm)                                                         | 6-10     |
| 2082-10  | 10" x 6" x 1.5" (25.4 cm x 15 cm x 3.8 cm)                                                         | 6-10     |
| 2081-21  | Additional tray insert with Scan-Guard™ silicone mat<br>10" x 6" x 1.5" (25.4 cm x 15 cm x 3.8 cm) | 10+      |

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible

Firm Sale – May be available for evaluation



SCANLAN®  
**V.I.P™ Sterilization Trays** (Valued Instrument Protection)  
Plastic



SCANLAN®  
**V.I.P™ Plastic Sterilization Trays**  
Scan-Guard™ Silicone Inserts, Detachable Lids

| Cat. No. | Description                                                                                                            | Capacity |
|----------|------------------------------------------------------------------------------------------------------------------------|----------|
| 2083-16  | Metal clasps, 13" x 7.5" x 0.7" (33.1 cm x 19.1 cm x 1.8 cm)                                                           | 10-15    |
| 2081-12  | Plastic clasps, 15" x 10" x 0.7" (38 cm x 25.4 cm x 1.8 cm)                                                            | 10-15    |
| 2084-12  | Metal clasps, additional tray insert and extra Scan-Guard™ silicone mat; 9.5" x 8" x 1.5" (24.1 cm x 20.3 cm x 3.8 cm) | 12-15    |
| 2082-20  | Plastic clasps, deep base<br>15" x 10" x 1.5" (38 cm x 25.4 cm x 3.8 cm)                                               | 10-15    |
| 2081-24  | Metal clasps, additional tray insert and Scan-Guard™ silicone mats<br>15" x 10" x 1.5" (38 cm x 25.4 cm x 3.8 cm)      | 20-30    |
| 2084-20  | Metal clasps, additional tray insert and Scan-Guard™ silicone mats<br>17.8" x 9.7" x 1.5" (45.2 cm x 25 cm x 3.8 cm)   | 24-30+   |

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible

Firm Sale – May be available for evaluation



SCANLAN®  
**Aorta Vein Punches**  
For Use in Myocardial Revascularization



Precisely engineered Scissor Cut™ blades featured on SCANLAN® single-use aorta/vein punches provide clean accurate incisions for myocardial revascularization. With large finger rings for greater control and accuracy, and tapered tips for an unobstructed view of the operating field, SCANLAN® single-use aorta/vein punches are available in two lengths (Original and Ultimate™), two styles (round and elliptical), and nine cutting sizes. Radiopaque. Not made with natural rubber latex.

SCANLAN®  
**Original Round**

| Cat. No. | Description |
|----------|-------------|
| 1001-600 | ● 2.7 mm    |
| 1001-601 | ● 3.5 mm    |
| 1001-602 | ● 4.0 mm    |
| 1001-604 | ● 4.5 mm    |
| 1001-605 | ● 5.0 mm    |
| 1001-606 | ● 5.5 mm    |
| 1001-607 | ● 6.0 mm    |

SCANLAN®  
**Ultimate Round**

| Cat. No. | Description |
|----------|-------------|
| 1001-620 | ● 2.7 mm    |
| 1001-621 | ● 3.5 mm    |
| 1001-622 | ● 4.0 mm    |
| 1001-624 | ● 4.5 mm    |
| 1001-626 | ● 5.0 mm    |
| 1001-628 | ● 5.5 mm    |
| 1001-629 | ● 6.0 mm    |

SCANLAN®  
**Original Elliptical**

| Cat. No. | Description    |
|----------|----------------|
| 1001-608 | ● 2.0 x 4.0 mm |
| 1001-610 | ● 2.0 x 6.0 mm |

SCANLAN®  
**Ultimate Elliptical**

| Cat. No. | Description    |
|----------|----------------|
| 1001-609 | ● 2.0 x 4.0 mm |
| 1001-611 | ● 2.0 x 6.0 mm |

1 per sterile package, 5 packages per box




SCANLAN®  
**A/C LOCATOR® & RADIOMARK®**  
Graft Markers



Scanlan offers two styles of single-use graft markers. U-shaped A/C LOCATOR® graft markers help to assess graft patency during postoperative cardiac angiography. RADIOMARK® graft markers are smaller, disk-shaped graft markers designed for use with peripheral and coronary grafts. Following graft placement, all graft markers can be applied using 5-0 or 6-0 sutures and will not migrate. SCANLAN® graft markers are flexible, lightweight, radiopaque and tissue-compatible to meet all implant criteria. Not made with natural rubber latex.

| Cat. No. | Description               |
|----------|---------------------------|
| 1001-83  | A/C LOCATOR® graft marker |
| 1001-831 | RADIOMARK® graft marker   |





# Vascular Tunneling System

## For Implementation of Subcutaneous Grafts

The SCANLAN® Vascular Tunneling System combines high quality grasping forceps with single-use tunneler sheaths for fast, accurate, atraumatic implantation of subcutaneous grafts.

SCANLAN®

### Grasping Forceps

The grasping forceps features include ...

- A choice of three styles: Vascular tunneler grasping forceps (20.5" / 52 cm), A/V Access tunneler grasping forceps (12" / 30.5 cm), and Full Curve A/V Access tunneler grasping forceps (12" / 30.5 cm)
- Each style will accommodate two sheath diameters (large or small)
- A knurled handle for positive control
- A handle indicator to verify tunneler curve direction
- Stainless steel construction, offering strength and long life, as well as quick, reliable, atraumatic graft placement
- Specially designed, stainless steel sterilization/storage case also available (see page 20)

SCANLAN®

### Grasping Forceps

Vascular  
Tunneler  
9009-24  
20.5" (52 cm)



A/V Access  
Tunneler  
9009-14  
12" (30.5 cm)



A/V Access  
Tunneler  
Partial Curve  
9009-13  
12" (30.5 cm)



A/V Access  
Tunneler  
Full Curve  
9009-15  
12" (30.5 cm)



# Vascular Tunneling System

## For Implementation of Subcutaneous Grafts

SCANLAN®

### Single-Use Sheaths

The single-use sheath features include ...

- Two sheath diameters (Large and Small) designed to fit each forceps style
- Allows placement of grafts without kinking, twisting, stretching, or pulling
- Smooth polymer sheaths and bullet tips offer near frictionless graft placement, essential for many new graft surfaces
- Minimal sheath diameters avoid excessive tissue disruption
- Color-coding provides easy length and diameter identification
- Radiopaque and latex-safe

SCANLAN®

### Single-Use Sheaths

1 per sterile package, 5 packages per box

9009-16  
20.5" (52 cm)



9009-18  
20.5" (52 cm)



9009-20  
12" (30.5 cm)



9009-22  
12" (30.5 cm)



9009-20  
12" (30.5 cm)



9009-22  
12" (30.5 cm)



9009-21  
12" (30.5 cm)



9009-23  
12" (30.5 cm)



Interior  
Diameters  
(mm)

For Grafts  
Up To  
(mm)

12

10

9.5

8

9

8

7

6

9

8

7

6

9

8

7

6




SCANLAN®

# Vascular Tunneling System

## For Implementation of Subcutaneous Grafts

Pictured  
9009-11 Set with Instrument tray 2150-03



SCANLAN®

### Vascular Tunneler and A/V Access Tunneler Grasping Forceps Set..... 9009-11

| Qty. | Cat. No. | Description                                                                     | Dimensions                                        |
|------|----------|---------------------------------------------------------------------------------|---------------------------------------------------|
| 1    | 9009-13  | Partial Curve A/V Access Tunneler                                               | 12"<br>(30.5 cm)                                  |
| 1    | 9009-14  | A/V Access Tunneler                                                             | 12"<br>(30.5 cm)                                  |
| 1    | 9009-15  | Full Curve A/V Access Tunneler                                                  | 12"<br>(30.5 cm)                                  |
| 1    | 9009-24  | Vascular Tunneler                                                               | 20.5"<br>(52 cm)                                  |
| 1    | 2150-03  | SCANLAN® V.I.P.™<br>Large capacity tray; Scan-Guard™<br>silicone brackets & mat | 29" x 8.4" x 1.5"<br>(73.7 cm x 21.3 cm x 3.8 cm) |

Set does not include tunneler sheaths  
All trays are steam, Immediate Use Steam Sterilization (IUSS), ETO, and plasma sterilization compatible  
*Firm Sale – May be available for evaluation*



SCANLAN®

# Vascular Tunneling System

## For Implementation of Subcutaneous Grafts

Pictured  
9009-10 Set with Instrument tray 9009-26



SCANLAN®

### Vascular Tunneler and A/V Access Tunneler Grasping Forceps Set..... 9009-10

| Qty. | Cat. No. | Description                                      | Dimensions                                     |
|------|----------|--------------------------------------------------|------------------------------------------------|
| 1    | 9009-14  | A/V Access Tunneler                              | 12"<br>(30.5 cm)                               |
| 1    | 9009-24  | Vascular Tunneler                                | 20.5"<br>(52 cm)                               |
| 1    | 9009-26  | SCANLAN® V.I.P.™<br>Storage / Sterilization Tray | 3.5" x 26" x 1.5"<br>(8.9 cm x 66 cm x 3.8 cm) |

Set does not include tunneler sheaths

SCANLAN®

### Vascular Tunneler and A/V Access Tunneler Grasping Forceps Set..... 9009-12

| Qty. | Cat. No. | Description                                      | Dimensions                                     |
|------|----------|--------------------------------------------------|------------------------------------------------|
| 1    | 9009-24  | Vascular Tunneler                                | 20.5"<br>(52 cm)                               |
| 1    | 9009-26  | SCANLAN® V.I.P.™<br>Storage / Sterilization Tray | 3.5" x 26" x 1.5"<br>(8.9 cm x 66 cm x 3.8 cm) |

Set does not include tunneler sheaths

All trays are steam, Immediate Use Steam Sterilization (IUSS), ETO, and plasma sterilization compatible  
*Firm Sale – May be available for evaluation*






# SCANLAN® VASCU-STATT® Single-Use Bulldog Clamps

Vascu-Stat® single-use bulldog clamps are a lightweight, cost-effective alternative to metal bulldog clamps.

Unlike metal bulldog clamps (clamping pressure weakens after repeated use), Vascu-Stat® single-use bulldog clamps offer a preset tension – ensuring consistent clamping pressure with each use.

The entire line of Vascu-Stat® single-use bulldog clamps feature...

- **Two designs** – straight or angled
- **Lightweight, cost-effective, plastic construction**
- **A wide range** of preset, consistent clamping pressures
- **Color-coded** – by clamp size and pressure
- **Single-use, sterile packaged**
- **Atraumatic** – jaw designs
- **Versatile** – may be applied by hand or with the specially designed VASCU-STATT® appliers
- **Radiopaque**
- **Not made with natural rubber latex**



## SCANLAN® Vascu-Stat® System

**Vascu-Stat®**  
Original bulldog clamp design

**Vascu-Stat® II**  
Firmier clamping pressure than original Vascu-Stat®

**Vascu-Stat® Plus**  
Vascu-Stat® PLUS single-use bulldog clamps offer the same benefits as the Vascu-Stat® original clamps, but features a cushioned Gentle Jaw™ surface that absorbs and softens as it provides an atraumatic occlusion of the vessel.

**Vascu-Stat® Micro**  
Micro Vascu-Stat® clamps offer even smaller clamp sizes – preferred for use in peripheral micro vascular and neuro surgery, as well as in many types of plastic and hand surgery.

Micro Vascu-Stat® single-use bulldog clamps also offer a radiopaque tag line for safe and easy visualization within the operating field.

## SCANLAN® Micro Vascu-Stat®


One size, available with straight jaws.  
Color-coded Purple.



**Micro Vascu-Stat®**  
1001-499 Micro, straight (purple)

## SCANLAN® Micro Vascu-Stat® II

One size, available with straight jaws.  
Color-coded Black.



**Micro Vascu-Stat® II**  
1001-529 Micro, straight (black)

2 each per sterile package, 5 packages per box



# SCANLAN® VASCU-STATT® Single-Use Bulldog Clamps

## SCANLAN® Vascu-Stat® & Vascu-Stat® Plus

Three sizes, available with straight or angled jaws.  
Color-coded.

|                               |  |                         |                        |
|-------------------------------|--|-------------------------|------------------------|
| Clamping Pressure 10-15 grams |  | <b>Vascu-Stat®</b>      |                        |
|                               |  | 1001-502                | Mini, straight (red)   |
|                               |  | 1001-503                | Mini, angled 45° (red) |
|                               |  | <b>Vascu-Stat® Plus</b> |                        |
|                               |  | 1001-552                | Mini, straight (red)   |
|                               |  | 1001-553                | Mini, angled 45° (red) |

|                               |  |                         |                          |
|-------------------------------|--|-------------------------|--------------------------|
| Clamping Pressure 25-30 grams |  | <b>Vascu-Stat®</b>      |                          |
|                               |  | 1001-500                | Midi, straight (white)   |
|                               |  | 1001-501                | Midi, angled 45° (white) |
|                               |  | <b>Vascu-Stat® Plus</b> |                          |
|                               |  | 1001-550                | Midi, straight (white)   |
|                               |  | 1001-551                | Midi, angled 45° (white) |

|                               |  |                         |                         |
|-------------------------------|--|-------------------------|-------------------------|
| Clamping Pressure 50-60 grams |  | <b>Vascu-Stat®</b>      |                         |
|                               |  | 1001-504                | Maxi, straight (blue)   |
|                               |  | 1001-505                | Maxi, angled 45° (blue) |
|                               |  | <b>Vascu-Stat® Plus</b> |                         |
|                               |  | 1001-554                | Maxi, straight (blue)   |
|                               |  | 1001-555                | Maxi, angled 45° (blue) |

2 per sterile package, 5 packages per box



## SCANLAN® Vascu-Stat® II & Vascu-Stat® II Plus

Three sizes, available with straight or angled jaws.  
Color-coded.

|                               |  |                            |                          |
|-------------------------------|--|----------------------------|--------------------------|
| Clamping Pressure 20-25 grams |  | <b>Vascu-Stat® II</b>      |                          |
|                               |  | 1001-530                   | Mini, straight (green)   |
|                               |  | 1001-531                   | Mini, angled 45° (green) |
|                               |  | <b>Vascu-Stat® II Plus</b> |                          |
|                               |  | 1001-570                   | Mini, straight (green)   |
|                               |  | 1001-571                   | Mini, angled 45° (green) |

|                               |  |                            |                           |
|-------------------------------|--|----------------------------|---------------------------|
| Clamping Pressure 75-80 grams |  | <b>Vascu-Stat® II</b>      |                           |
|                               |  | 1001-532                   | Midi, straight (yellow)   |
|                               |  | 1001-533                   | Midi, angled 45° (yellow) |
|                               |  | <b>Vascu-Stat® II Plus</b> |                           |
|                               |  | 1001-572                   | Midi, straight (yellow)   |
|                               |  | 1001-573                   | Midi, angled 45° (yellow) |

|                                 |  |                            |                           |
|---------------------------------|--|----------------------------|---------------------------|
| Clamping Pressure 165-175 grams |  | <b>Vascu-Stat® II</b>      |                           |
|                                 |  | 1001-534                   | Maxi, straight (orange)   |
|                                 |  | 1001-535                   | Maxi, angled 45° (orange) |
|                                 |  | <b>Vascu-Stat® II Plus</b> |                           |
|                                 |  | 1001-574                   | Maxi, straight (orange)   |
|                                 |  | 1001-575                   | Maxi, angled 45° (orange) |

SCANLAN®  
**VASCU-STATT®**  
Single-Use Bulldog Clamps

SCANLAN®

**Vascu-Stat® Approximator & Vascu-Stat® Plus Approximator**

Used to facilitate end-to-end anastomosis. Three sizes, available with straight or angled jaws. Color-coded.

Clamping Pressure 10-15 grams



**Vascu-Stat® Approximator**  
1001-507 Mini, straight (red)  
1001-508 Mini, angled 45° (red)  
**Vascu-Stat® Plus Approximator**  
1001-560 Mini, straight (red)  
1001-561 Mini, angled 45° (red)

Not Shown Actual Size


Clamping Pressure 25-30 grams



**Vascu-Stat® Approximator**  
1001-509 Midi, straight (white)  
1001-510 Midi, angled 45° (white)  
**Vascu-Stat® Plus Approximator**  
1001-562 Midi, straight (white)  
1001-563 Midi, angled 45° (white)

Not Shown Actual Size

Clamping Pressure 50-60 grams



**Vascu-Stat® Approximator**  
1001-511 Maxi, straight (blue)  
1001-512 Maxi, angled 45° (blue)  
**Vascu-Stat® Plus Approximator**  
1001-564 Maxi, straight (blue)  
1001-565 Maxi, angled 45° (blue)

Not Shown Actual Size

1 per sterile package, 5 packages per box



SCANLAN®

**Vascu-Stat® II Approximator & Vascu-Stat® II Plus Approximator**

Used to facilitate end-to-end anastomosis. Three sizes, available with straight or angled jaws. Color-coded.


Clamping Pressure 20-25 grams



**Vascu-Stat® II Approximator**  
1001-540 Mini, straight (green)  
1001-541 Mini, angled 45° (green)  
**Vascu-Stat® II Plus Approximator**  
1001-580 Mini, straight (green)  
1001-581 Mini, angled 45° (green)

Not Shown Actual Size

Clamping Pressure 75-80 grams



**Vascu-Stat® II Approximator**  
1001-542 Midi, straight (yellow)  
1001-543 Midi, angled 45° (yellow)  
**Vascu-Stat® II Plus Approximator**  
1001-582 Midi, straight (yellow)  
1001-583 Midi, angled 45° (yellow)

Not Shown Actual Size

Clamping Pressure 165-175 grams



**Vascu-Stat® II Approximator**  
1001-544 Maxi, straight (orange)  
1001-545 Maxi, angled 45° (orange)  
**Vascu-Stat® II Plus Approximator**  
1001-584 Maxi, straight (orange)  
1001-585 Maxi, angled 45° (orange)

Not Shown Actual Size

SCANLAN®

**VASCU-STATT®**  
Single-Use Bulldog Clamps

SCANLAN®

**Vascu-Stat® Appliers**

Vascu-Stat® appliers are designed to be used exclusively with Vascu-Stat® single-use bulldog clamps. They are a convenient, easy way to apply Vascu-Stat® single-use bulldog clamps — especially in deep, hard to reach surgical fields.

| Cat. No.  | Description                                            |
|-----------|--------------------------------------------------------|
| 1001-506  | VASCU-STATT® Applier, Straight, 8 1/2" (21.5 cm)       |
| 1001-506A | VASCU-STATT® Applier, Angled, 8 1/2" (21.5 cm)         |
| 1001-506M | Micro VASCU-STATT® Applier, Straight, 8 1/2" (21.5 cm) |
| 9009-506  | VASCU-STATT® Applier, Straight, 15 1/2" (39.4 cm)      |
| 9009-506A | VASCU-STATT® Applier, Angled, 15 1/2" (39.4 cm)        |



1001-506  
VASCU-STATT®  
Applier, Straight



1001-506A  
VASCU-STATT®  
Applier, Curved



1001-506M  
Micro VASCU-STATT®  
Applier, Straight



9009-506  
CardioVasive™  
VASCU-STATT®  
Applier, Straight



9009-506A  
CardioVasive™  
VASCU-STATT®  
Applier, Curved





SCANLAN®  
**SURG-I-LOOP®**  
Vessel Loops

Designed to assist surgeons in complicated surgical fields by providing retraction, occlusion and identification of arteries, veins, nerves and ureters.

- **Unique design** – provides access to complicated surgical fields
- **Versatile** – can retract and occlude arteries, veins, nerves and ureters
- **Durable** – retains preset tension
- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description      |
|----------|------------------|
| 1001-75  | Red, Mini        |
| 1001-76  | Red, Maxi        |
| 1001-77  | Blue, Mini       |
| 1001-78  | Blue, Maxi       |
| 1001-78A | Blue, Super Maxi |
| 1001-79  | Yellow, Mini     |
| 1001-80  | Yellow, Maxi     |
| 1001-81  | White, Mini      |
| 1001-81M | White, Micro     |
| 1001-82  | White, Maxi      |

Length 16 1/2" (42 cm)

2 per sterile package,  
10 packages per box


SCANLAN®  
**SURG-I-LOOP® PLUS**

SURG-I-LOOP® *Plus* features a solid stainless steel tip for fast, accurate placement / passing in complicated surgical fields. Ideal for minimally-invasive cardiac surgery. Radiopaque. Not made with natural rubber latex.

| Cat. No. | Description                                  |
|----------|----------------------------------------------|
| 1001-81P | Blunt needle (1.2 mm)<br>(light green)       |
| 1001-82P | Semi-sharp needle (1.2 mm)<br>(light orange) |

2 per sterile package,  
10 packages per box

|   |                          |
|---|--------------------------|
| ● | Super Maxi (5 x 1.27 mm) |
| ● | Maxi (2.5 x 1.02 mm)     |
| ● | Mini (1.5 x .88 mm)      |
| ● | Micro (0.5 x .26 mm)     |



SCANLAN®  
**CLAMP JAW COVERS**



SCANLAN®  
**SOFT-GRIP™**  
Fabric Clamp Covers

Special foam and fabric combination provides a soft, non-slip clamp surface.

- **Special material** – assures no fibrous residue in operating field
- **Versatile** – designed to fit a variety of instruments
- **Ease of application** – designed with convenient pull tabs
- **Adaptable** – may be trimmed to desired length
- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description                                                               |
|----------|---------------------------------------------------------------------------|
| 1001-740 | SOFT-GRIP™ fabric clamp covers<br>Length 5" (127 mm), Width 3/8" (9.5 mm) |

2 per sterile package,  
10 packages per box

SCANLAN®  
**SURG-I-PAW®**  
Instrument Jaw Covers

Designed for Intestinal, Thoracic, Pediatric, and other specialty types of surgery.

- **Atraumatic** – converts standard clamps into atraumatic instruments
- **Deep serrations** – provides firm, atraumatic occlusion
- **Contoured design** – provides secure pressure – will not slip or twist
- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description           |
|----------|-----------------------|
| 1001-72  | 1/10" (2.5 mm), Red   |
| 1001-73  | 2/10" (5.0 mm), White |
| 1001-74  | 3/10" (7.5 mm), Blue  |

2 per sterile package,  
10 packages per box



SCANLAN®  
SINGLE-USE



SCANLAN®  
**PINCAP®**  
Orthopedic Pin Protectors

Unique single-use protective covers for exposed tips of orthopedic pins and wires.

- **Unique design** – provides excellent protection from dermal abrasions, damage to bedding, clothing snags, and surgical glove snags
- **Size selection** – designed to fit most common pins and wires
- **Double-wall design** – provides additional protection
- **Color coded** – providing easy identification
- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description                                                                               |
|----------|-------------------------------------------------------------------------------------------|
| 301-0110 | White (fits 3 sizes K-wire)<br>0.028" (0.71 mm),<br>0.035" (0.89 mm),<br>0.045" (1.14 mm) |
| 301-0120 | Pink and White<br>0.062" (1.57 mm)                                                        |
| 301-0130 | Blue and White<br>0.078" (5/64") (1.98 mm)                                                |

2 per sterile package,  
20 packages per box



SCANLAN®  
**SUTURE BOOT™**  
Instrument Jaw Covers

Designed to facilitate the grasping and tagging of suture.

- **Atraumatic** – firmly holds the finest suture without damage
- **Secure** – will not slip
- **Color coded** – bright yellow and pink for immediate identification
- **Foam cartridge** – designed for fast / easy access
- **Radiopaque**
- **Not made with natural rubber latex**



| Cat. No. | Description                                         |
|----------|-----------------------------------------------------|
| 1001-513 | Original Mini, Yellow<br>(3 pair per cartridge)     |
| 1001-514 | Original Standard, Yellow<br>(5 pair per cartridge) |
| 1001-515 | Large Mini, Pink<br>(3 pair per cartridge)          |
| 1001-516 | Large Standard, Pink<br>(5 pair per cartridge)      |

1 per sterile package,  
5 packages per box



SCANLAN®  
SINGLE-USE



SCANLAN®  
**SUTURE-STATT®**  
Suture Tags

Designed to facilitate the grasping and tagging of suture during ophthalmic or delicate procedures.

- **Secure** – firmly holds the finest suture without damage
- **Keeps sutures out of the operative field**
- **Lightweight** – only 0.5 grams
- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description |
|----------|-------------|
| 1001-520 | Suture Tags |

2 per sterile package,  
5 packages per box

SCANLAN®  
**SOLEM® Vein Holder**

Ideal for assisting surgeons during vessel anastomosis, the SCANLAN® Solem™ vein holder gently yet firmly grasps vein grafts and precisely positions suture without trauma. Its slim profile and concave tip assures optimal exposure and the lightly serrated clip provides steady, consistent pressure, minimizing graft movement.

- **Radiopaque**
- **Not made with natural rubber latex**

| Cat. No. | Description                          |
|----------|--------------------------------------|
| 1001-760 | Solem™ vein holder, 6 1/2" (16.5 cm) |

1 per sterile package,  
5 packages per box

SCANLAN®  
SURGICAL SKIN MARKERS

SCANLAN®  
**PrepMark™**

Available in gentian violet with regular tip, SCANLAN® PrepMARK™ surgical skin markers are made with a unique prep-resistant ink and ideal for marking, identifying, coding, and contouring the surgical area. Latex-safe.

Regular  
Tip

| Cat. No. | Description                                                                                                     |
|----------|-----------------------------------------------------------------------------------------------------------------|
| 1001-730 | PrepMARK™ Surgical Skin Marker, Sterile, plus ruler<br>Box of 50 (1 per sterile package with ruler, 50 per box) |
| 1001-731 | PrepMARK™ Surgical Skin Marker, Non-Sterile<br>Container of 50                                                  |



Available in gentian violet with both regular and extra fine tips, SCANLAN® surgical skin markers are appropriate for marking, identifying, coding, and contouring the surgical area. Not made with natural rubber latex.

| Regular Tip | Extra Fine Tip | Description    |
|-------------|----------------|----------------|
| Cat. No.    | Cat. No.       |                |
| 1001-710    | 1001-725       | Gentian Violet |


1 per sterile package with ruler,  
5 packages per box






# SINGLE-USE & INSTRUMENT CARE Products

INCHES



CENTIMETERS



For more information  
or to place an order,  
please contact your local Scanlan  
representative,  
or call...

**U.S. & Canada:** 800-328-9458

**International:** 651-298-0997

**FAX:** 651-298-0018

[www.scanlaninternational.com](http://www.scanlaninternational.com)

