


SCANLAN®

Temporary Occlusion Products

VASCU-STATT®

Single-Use Bulldog Clamps



RELIANCE BULLDOG CLAMPS

Reliable Surgeon-Controlled Vessel Clamps




YASARGIL

Reuseable Temporary Occlusion Clips



HEIFETZ™

Reuseable Temporary Occlusion Clips



100 Years of Surgical Excellence... and Our Journey Continues



Joseph
Scanlan

Russell
Scanlan

Tim
Scanlan

Adam
Scanlan

Chenoa
Scanlan

Since 1921, The Scanlan family has supplied hospital operating rooms with the highest quality surgical instrumentation and products. We are dedicated to work hand-in-hand with surgeons and perioperative nurses around the world to develop new and innovative products.

We offer a wide selection of instruments designed specifically for temporary occlusion – Each reflecting Scanlan’s commitment to quality through superior designs, the highest quality medical-grade materials, rigorous evaluation and testing.

Our goal is simple – provide the finest surgical products to assist the world’s surgeons, nurses and ultimately... their patients.

Sincerely,

Timothy M. Scanlan
President and Chief Executive Officer

Russell T. Scanlan
Chief Operating Officer



VASCU-STAT®

Single-Use Bulldog Clamps

VASCU-STAT® single-use bulldog clamps are a lightweight, cost-effective alternative to metal bulldog clamps.

Features Include:

- Angled and straight designs
- Variety of clamping pressures
- Color coded for easy identification
- Atraumatic jaws
- Apply by hand or with Scanlan Vascu-Stat® appliers
- Radiopaque
- Not made with natural rubber latex

VASCU-STAT® bulldog clamps feature a consistent pre-set tension unlike metal bulldog clamps where pressure can weaken after repeated use.



Parallel clamp designed for
maximum occlusion and
atraumatic pressure



SCANLAN®
VASCU-STATT®
Single-Use Bulldog Clamps

SCANLAN®
VASCU-STATT®
Regular Clamping Pressure

Original bulldog clamp design.



Cat. No.	Description
1001-502	Mini, Straight (Red) 
1001-503	Mini, Angled 45° (Red) 
1001-500	Midi, Straight (White) 
1001-501	Midi, Angled 45° (White) 
1001-504	Maxi, Straight (Blue) 
1001-505	Maxi, Angled 45° (Blue) 

2 per sterile package, 5 packages per box.



SCANLAN®
VASCU-STATT® Approximator
Regular Clamping Pressure

Facilitates end-to-end anastomosis.

Cat. No.	Description
1001-507	Mini, Straight (Red) 
1001-508	Mini, Angled 45° (Red) 
1001-509	Midi, Straight (White) 
1001-510	Midi, Angled 45° (White) 
1001-511	Maxi, Straight (Blue) 
1001-512	Maxi, Angled 45° (Blue) 

SCANLAN®
VASCU-STATT® Plus
Single-Use Bulldog Clamps

SCANLAN®
VASCU-STATT® Plus
Regular Clamping Pressure

VASCU-STATT® PLUS offers the same benefits as VASCU-STATT® clamps, but feature a cushioned Gentle Jaw™ surface that absorbs and softens as it provides an atraumatic occlusion of the vessel.

Cat. No.	Description
1001-552	Mini, Straight Gentle Jaw™ (Red) 
1001-553	Mini, Angled 45° Gentle Jaw™ (Red) 
1001-550	Midi, Straight Gentle Jaw™ (White) 
1001-551	Midi, Angled 45° Gentle Jaw™ (White) 
1001-554	Maxi, Straight Gentle Jaw™ (Blue) 
1001-555	Maxi, Straight Gentle Jaw™ (Blue) 

2 per sterile package, 5 packages per box.



SCANLAN®
VASCU-STATT® Plus Approximator
Regular Clamping Pressure

Facilitates end-to-end anastomosis.

Cat. No.	Description
1001-560	Mini, Straight Gentle Jaw™ (Red) 
1001-561	Mini, Angled 45° Gentle Jaw™ (Red) 
1001-562	Midi, Straight Gentle Jaw™ (White) 
1001-563	Midi, Angled 45° Gentle Jaw™ (White) 
1001-564	Maxi, Straight Gentle Jaw™ (Blue) 
1001-565	Maxi, Straight Gentle Jaw™ (Blue) 

SCANLAN®
VASCU-STATT® II
Single-Use Bulldog Clamps

SCANLAN®
VASCU-STATT® II
Firm Clamping Pressure

Original bulldog clamp design.

Cat. No.	Description	
1001-530	Mini, Straight (Green)	
1001-531	Mini, Angled 45° (Green)	
1001-532	Midi, Straight (Yellow)	
1001-533	Midi, Angled 45° (Yellow)	
1001-534	Maxi, Straight (Orange)	
1001-535	Maxi, Angled 45° (Orange)	

2 per sterile package, 5 packages per box.



SCANLAN®
VASCU-STATT® II Plus
Single-Use Bulldog Clamps

SCANLAN®
VASCU-STATT® II Plus Approximator
Firm Clamping Pressure

Facilitates end-to-end anastomosis.

SCANLAN®
VASCU-STATT® II Plus
Firm Clamping Pressure

VASCU-STATT® II PLUS offers the same benefits as VASCU-STATT® clamps, but feature a cushioned Gentle Jaw™ surface that absorbs and softens as it provides an atraumatic occlusion of the vessel.

Cat. No.	Description	
1001-570	Mini, Straight Gentle Jaw™ (Green)	
1001-571	Mini, Angled 45° Gentle Jaw™ (Green)	
1001-572	Midi, Straight, Gentle Jaw™ (Yellow)	
1001-573	Midi, Angled 45° Gentle Jaw™ (Yellow)	
1001-574	Maxi, Straight Gentle Jaw™ (Orange)	
1001-575	Maxi, Angled 45° Gentle Jaw™ (Orange)	

2 per sterile package, 5 packages per box.



SCANLAN®
VASCU-STATT® II Approximator
Firm Clamping Pressure

Facilitates end-to-end anastomosis.

Cat. No.	Description	
1001-540	Mini, Straight (Green)	
1001-541	Mini, Angled 45° (Green)	
1001-542	Midi, Straight (Yellow)	
1001-543	Midi, Angled 45° (Yellow)	
1001-544	Maxi, Straight (Orange)	
1001-545	Maxi, Angled 45° (Orange)	

SCANLAN®
MICRO VASCU-STATT®
Single-Use Bulldog Clamps

SCANLAN®
VASCU-STATT® Micro

Micro VASCU-STATT® clamps offer smaller clamp sizes – preferred for use in peripheral micro vascular and neuro surgery, as well as in many types of plastic and hand surgery. Micro VASCU-STATT® single-use bulldog clamps also offer a radiopaque tag line for safe and easy visualization within the operating field.

SCANLAN®
VASCU-STATT® Micro
Regular Clamping Pressure

Cat. No.	Description
1001-499	Straight (Purple) 

2 per sterile package, 5 packages per box.

SCANLAN®
VASCU-STATT® II Micro
Firm Clamping Pressure

Cat. No.	Description
1001-529	Straight (Black) 

SCANLAN®
VASCU-STATT®
Single-Use Bulldog Clamps

SCANLAN®
VASCU-STATT® Appliers

VASCU-STATT® appliers are designed to be used exclusively with VASCU-STATT® single-use bulldog clamps. They are a convenient, easy way to apply VASCU-STATT® single-use bulldog clamps — especially in deep, hard to reach surgical fields.



Cat. No.	Description	Total Length
1001-506	 VASCU-STATT® Applier Straight	8-1/2" (21.5 cm)
1001-506A	 VASCU-STATT® Applier, Angled	
1001-506M	 Micro VASCU-STATT® Applier, Straight	

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.



SCANLAN®

Reliance Bulldog Clamps

Reliable Surgeon-Controlled Vessel Clamps

The SCANLAN® Reliance Bulldog Clamp has been designed to provide secure, robust clamping pressure to assist in the temporary occlusion of vessels. Manufactured with the highest quality stainless steel, its 1 x 2 DeBakey jaws ensure a firm yet atraumatic clamping pressure. A unique suture port allows optional suture tagging to enhance surgical location and recovery.

Features Include:

- 1 x 2 DeBakey jaws
- 350 grams closing pressure (approximate)
- Suture port
- Highest quality stainless steel
- Reusable / durable
- Low profile design
- Radiopaque



Cat. No.	Description
3795-52	40 mm Straight DeBakey Jaws 2-3/4" (7 cm) 
3795-50	25 mm Straight DeBakey Jaws 2" (5 cm) 
3795-56	15 mm Straight DeBakey Jaws 1-5/8" (4 cm) 
3795-54	7 mm Straight DeBakey Jaws 1-1/4" (3.2 cm) 

Cat. No.	Description
3795-51	25 mm Curved DeBakey Jaws 2" (5 cm) 
3795-57	15 mm Straight DeBakey Jaws 1-5/8" (4 cm) 

All SCANLAN® Reliance Bulldog Clips are re-usable and steam, flash (immediate use), ETO and plasma sterilization compatible. Not intended for implantation.



SCANLAN®


Reliance Bulldog Clamps

Reliable Surgeon-Controlled Vessel Clamps

The SCANLAN® Reliance Bulldog LP Clamp is based on the same platform as the current Reliance Bulldog Clamps but features a lower clamping pressure. The Low Pressure bulldog is designed to fit through a variety of introducer ports, making this a great option for use in situations when the patient's anatomy (or chosen surgical technique) is best served with lower pressure clamping. With its new color to differentiate this from other clamps, the SCANLAN® Reliance Bulldog LP Clamp sets the "gold standard" for effective clamping.

Reliance Bulldog LP Features

- Lower Clamping Pressure – provides approximately 205 grams of clamping pressure in initial use
- Expands Clamping Options for certain patients, techniques and anatomy
- Multiple Application Methods – including Scanlan Applier / Introducer
- Optimized Dimensions – length/height fits through commonly used ports and can be used with multiple instruments in the same port
- 18 mm Jaw Length – versatile, workhorse size
- Visually Differentiated – gold color stands out from other Reliance Bulldog clamps
- DeBakey Style Jaws – provides firm, atraumatic occlusion



Cat. No.	Description
3795-44	18 mm Straight Light Pressure DeBakey Jaws 1-3/4" (4.4 cm) 

All SCANLAN® Reliance Bulldog Clips are re-usable and steam, flash (immediate use), ETO and plasma sterilization compatible. Not intended for implantation.

PLEASE NOTE:
Not available in the European Union or Australia.
All products are not available in all regions. Please contact Scanlan customer service for product availability in your region.

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.



SCANLAN®


Reliance Bulldog Clamps

Reliable Surgeon-Controlled Vessel Clamps

SCANLAN®

Reliance Bulldog Appliers

- Custom jaw allows for secure deployment of bulldog into the surgical field
- Designed for use with 12 mm cannula



Cat. No.	Description		Working Length	Total Length
9009-393		Straight Jaws Shaft 8 mm With Ratchet	5" (12.7 cm)	9" (23 cm)
9009-395			13-3/4" (35 cm)	17-1/3" (44 cm)

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.



SCANLAN®

Reliance Bulldog Clamps

Reliable Surgeon-Controlled Vessel Clamps



SCANLAN®

Recommended Reliance Bulldog Starter Pack

Recommended	QTY	Cat. No.	Description	Total Length
SCANLAN® V.I.P.™ Sterilization Tray	Choose 1	2133-03	SCANLAN® V.I.P.™ Micro Sterilization Tray With Scan-Guard™ Silicone Mat - Includes Hardware (Holds 3-4 Clamps)	1-1/2" x 3-1/2" x 1" (3.81 cm x 8.89 cm x 2.54 cm)
		2135-03	SCANLAN® V.I.P.™ Reliance Bulldog Sterilization Tray With Scan-Guard™ Silicone Mat (Holds 8-10 Clamps)	8-1/4" x 3-1/4" x 1-1/2" (21 cm x 8.3 cm x 3.8 cm)
Applier/Introducer	Choose 1	9009-393	Straight Jaws With Ratchet	9" (23 cm)
		9009-395		17-1/3" (44.5 cm)
SCANLAN® Reliance Bulldog Clamps	Up To 4 Of Any Combination	3795-52	40 mm Straight DeBakey Jaws	2-3/4" (7 cm)
		3795-50	25 mm Straight DeBakey Jaws	2" (5 cm)
		3795-56	15 mm Straight DeBakey Jaws	1-5/8" (4 cm)
		3795-54	7 mm StraightDeBakey Jaws	1-1/4" (3.2 cm)
		3795-51	25 mm Curved DeBakey Jaws	2" (5 cm)
		3795-57	15 mm Curved DeBakey Jaws	1-5/8" (4 cm)



SCANLAN® Yasargil temporary occlusion clips offer surgeons a SMALL alternative to large, bulky metal or disposable clips with similar closing pressure ranges. A wide selection of clip designs provided to match specific surgical applications, such as internal carotid artery, internal mammary artery, radial arteries, etc.

SCANLAN® Yasargil temporary occlusion clips feature:

- Reusable / Durable – Allowing multiple sterilization cycles, resulting in low cost per procedure
- Low Profile Design – Providing excellent visualization
- Two Color-Coded Ranges of Closing Pressures – Offering reliable pre-set tensions
- Specially Designed Appliers – Providing controlled placement of clips
- Small Clip Applier Surface – Reducing suture entanglements
- Cobalt chromium alloy- High strength, resilient



SCANLAN®
Yasargil Clips
Light Pressure

Cat. No.	Description	
7 mm 3785-20	Slight Curve	
7 mm 3785-21	Straight	
9 mm 3785-22	Straight	
9 mm 3785-24	Slight Curve	
11 mm 3785-26	Straight	
11 mm 3785-27	Slight Curve	
15 mm 3785-41	Straight	
15 mm 3785-42	Slight Curve	
20 mm 3785-43	Straight	

SCANLAN®
Yasargil Clips
Firm Pressure

Cat. No.	Description	
7 mm 3786-20	Slight Curve	
7 mm 3786-21	Straight	
9 mm 3786-22	Straight	
9 mm 3786-23	Angled	
9 mm 3786-24	Slight Curve	
9 mm 3786-25	Strong Curve	

Cat. No.	Description	
11 mm 3786-26	Straight	
11 mm 3786-27	Slight Curve	
11 mm 3786-29	Strong Curve	
15 mm 3786-41	Straight	
15 mm 3786-42	Slight Curve	
20 mm 3786-43	Straight	

All SCANLAN® Yasargil Temporary Occlusion Clips are re-usable and steam, flash (immediate use), ETO and plasma sterilization compatible. Not intended for implantation.



SCANLAN®
Yasargil
Applying Applying Forceps – Standard

Cat. No.	Description	Total Length
3786-28	 Bayonet Style Flat Handle	7-1/8" (18.1 cm)



SCANLAN®

Yasargil Mini

Reuseable Temporary Occlusion Clips

SCANLAN® Yasargil temporary occlusion clips offer surgeons a SMALL alternative to large, bulky metal or disposable clips with similar closing pressure ranges. A wide selection of clip designs provided to match specific surgical applications, such as internal carotid artery, internal mammary artery, radial arteries, etc.

SCANLAN® Yasargil temporary occlusion clips feature:

- Reusable / Durable – Allowing multiple sterilization cycles, resulting in low cost per procedure
- Low Profile Design – Providing excellent visualization
- Specially Designed Appliers – Providing controlled placement of clips
- Small Clip Applier Surface – Reducing suture entanglements
- Cobalt chromium alloy- High strength, resilient



SCANLAN®

Yasargil Clips

Mini

	Cat. No.	Description	
4 mm	3786-30	Straight	
6 mm	3786-31	Straight	
7 mm	3786-32	Slight Curve	

All SCANLAN® Yasargil Temporary Occlusion Clips are re-usable and steam, flash (immediate use), ETO and plasma sterilization compatible. Not intended for implantation.



SCANLAN®

Yasargil Mini

Reuseable Temporary Occlusion Clips

3786-40



SCANLAN®

Yasargil

Applying Applying Forceps – Mini

Cat. No.	Description	Total Length
3786-40	 Bayonet Style Flat Handle	8-1/2" (21.5 cm)

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.




Reuseable Temporary Occlusion Clips

SCANLAN® Heifetz™ temporary occlusion clips offer surgeons a SMALL alternative to large, bulky metal or disposable clips with similar closing pressure ranges. A wide selection of clip designs provided to match specific surgical applications, such as internal carotid artery, internal mammary artery, radial arteries, etc.

SCANLAN® Heifetz™ temporary occlusion clips feature:

- Reusable / Durable – Allowing multiple sterilization cycles, resulting in low cost per procedure
- Low Profile Design – Providing excellent visualization
- Two Color-Coded Ranges of Closing Pressures – Offering reliable pre-set tensions
- Specially Designed Appliers – Providing controlled placement of clips
- Small Clip Applier Surface – Reducing suture entanglements
- Cobalt chromium alloy- High strength, resilient



SCANLAN®
Heifetz™ Clips
Light Pressure

Cat. No.	Description	
3795-04	Straight 8 mm Jaw Length 1-3/4" Jaw Width	
3795-06	Angled On Edge 8 mm Jaw Length 1-3/4" Jaw Width	
3795-08	Curved 8 mm Jaw Length 1-3/4" Jaw Width	
3795-20	Straight 12 mm Jaw Length 2-1/4" Jaw Width	
3795-22	Curved 12 mm Jaw Length 2-1/4" Jaw Width	


All SCANLAN® Heifetz™ Temporary Occlusion Clips are re-usable and steam, flash (immediate use), ETO and plasma sterilization compatible. Not intended for implantation.

SCANLAN®
Heifetz™ Clips
Firm Pressure

Cat. No.	Description	
3795-14	Straight 8 mm Jaw Length 1-3/4" Jaw Width	
3795-16	Angled On Edge 8 mm Jaw Length 1-3/4" Jaw Width	
3795-18	Curved 8 mm Jaw Length 1-3/4" Jaw Width	
3795-24	Straight 12 mm Jaw Length 2-1/4" Jaw Width	
3795-28	Curved 12 mm Jaw Length 2-1/4" Jaw Width	



Reuseable Temporary Occlusion Clips



3795-01

SCANLAN®
Heifetz™
Appliers

Cat. No.	Description	Total Length
3795-01	Straight	7-1/4" (18.5 cm)
3795-01A	Angled	

Please Note...
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.



SCANLAN®

Temporary Occlusion Products

For more information or to place an order,
please contact your local Scanlan
representative,
or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

www.scanlaninternational.com



MK-B-TOP Rev. B

SCANLAN INTERNATIONAL, Inc.
One Scanlan Plaza • Saint Paul, Minnesota 55107 • U.S.A. • 800-328-9458 • 651-298-0997 • FAX 651-298-0018
©2023 Scanlan International, Inc.