

SCANLAN® Microsurgical Instruments



Feel the Difference of Precision Craftsmanship



Highest Quality Surgical Instruments
Since 1921



Joseph
Scanlan

Adam
Scanlan

Tim
Scanlan

Russell
Scanlan

Chenoa
Scanlan

Ron
Peick

Dear Customer,

Since 1921, Scanlan International, Inc. has been designing and manufacturing the highest quality surgical instrumentation. Our Microsurgical line is a unique collection of instruments for procedures working within small to medium sized surgical fields. Each instrument is designed with the highest grade stainless steel or titanium alloys, and are handcrafted with meticulous precision... allowing excellent tactile response and control.

SCANLAN® Super-Micro Fine instrumentation is designed to meet the needs of all microsurgeons who work under a high power operating microscope or utilize surgical loupes to facilitate microvascular surgical technique by anastomosing small nerves and vessels.

Sincerely,

Timothy M. Scanlan
President and Chief Executive Officer

SCANLAN®

Super-Micro Fine Instrumentation

Highest Quality Materials and Manufacturing to Provide a Difference You Can Feel.

SCANLAN® Super-Micro Fine Needle Holders

- High Quality Titanium and Stainless Steel
- Titanium Needle Holders are Non-magnetic – eliminates needle jump
- Titanium Needle Holders are Lightweight – reduces the chance for wrist fatigue
- Low Profile Tapered Shanks for improved visibility under the microscope
- Designed For Use With 9-0 Suture (100 micron needles) or smaller
- 0.3 mm **Super-Micro Fine** tip
- **DIAMOND DUST™** jaws provides a secure grip while protecting the needle from burring
- Non-locking



SCANLAN® Super-Micro Fine Forceps

- High Quality Titanium and Stainless Steel
- Titanium Forceps are Non-magnetic – eliminates needle from clinging to instrument
- Titanium Forceps are Lightweight – reduces the chance for wrist fatigue
- Low Profile Tapered Shanks for improved visibility under the microscope
- Specifically designed for manipulation of small vessels
- 0.3 mm **Super-Micro Fine** Tip
- **DIAMOND DUST™** coating provides superior grasping of tissue



SCANLAN® Super-Micro Fine Scissors

- High Quality Hardened Stainless Steel for industry leading durability
- **SCANLAN® Premier™** scissors use proprietary technology that enhances the hardness and corrosion resistant properties of stainless steel
- 100% Cut Tested to ensure a high level of quality
- 12-month Guarantee on the sharpness of Scanlan® **Premier™** scissors



All Scanlan Instruments are backed by a 10-year Guarantee against manufacturing material defects



SCANLAN®

MICROSURGICAL NEEDLE HOLDERS

SCANLAN®
Super-Micro Fine Needle Holders
for 9-0 Suture

3303-646



3303-622



Cat. No.	Description	Handle	Recommended Suture Size	Total Length
NEW 3303-646	SCANLAN® <i>Super-Micro Fine</i> Titanium Needle Holder Curved, 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Round	9-0 and Smaller	5" (12.7 cm)
3303-603	SCANLAN® <i>Super-Micro Fine</i> Titanium Needle Holder Straight, 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Flat		
NEW 3303-664	SCANLAN® <i>Super-Micro Fine</i> Titanium Needle Holder Straight, 0.3 mm DIAMOND DUST™ Jaws, Regular Box Lock Without Ratchet	Round		7 1/4" (18.5 cm)
NEW 3303-622	SCANLAN® <i>Super-Micro Fine</i> Titanium Needle Holder Straight, 0.3 mm DIAMOND DUST™ Jaws, Regular Box Lock Without Ratchet	Flat		
3303-667	SCANLAN® <i>Super-Micro Fine</i> Titanium Needle Holder Curved, 0.3 mm DIAMOND DUST™ Jaws, Regular Box Lock Without Ratchet	Round		



SCANLAN®

MICROSURGICAL NEEDLE HOLDERS

SCANLAN®
Micro Needle Holders
for 8-0 Suture

3003-177-4



Cat. No.	Description	Handle	Recommended Suture Size	Total Length
3003-177-3	SCANLAN® Titanium Jacobson Micro Needle Holder Curved, 0.3 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	Round	8-0 and Smaller	7 1/4" (18.5 cm)
3003-177-4	SCANLAN® Titanium Jacobson Micro Needle Holder Curved, 0.4 mm DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet			

SCANLAN®
Micro Needle Holders
for 7-0 Suture

3003-207



Cat. No.	Description	Handle	Recommended Suture Size	Total Length
6006-44	SCANLAN® Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	Flat	7-0 and Smaller	4 3/4" (12 cm)
3003-207	SCANLAN® Titanium Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet			5 3/4" (14.5 cm)




SCANLAN®

MICROSURGICAL NEEDLE HOLDERS

SCANLAN®

Micro Needle Holders

for 6-0 Suture




Cat. No.	Description	Handle	Recommended Suture Size	Total Length
3003-452	SCANLAN® Titanium Richards Micro Needle Holder Extra Fine Straight DIAMOND DUST™ Jaws With Ratchet	Ring	6-0 and Smaller	6 3/4" (17 cm)

SCANLAN®

Micro Needle Holders

for 5-0 Suture



Cat. No.	Description	Handle	Recommended Suture Size	Total Length
4155-07	SCANLAN® Webster Needle Holder Straight delicate jaws Offset bows	Ring	5-0 and Smaller	5 1/8" (13 cm)
4155-06	SCANLAN® Webster Needle Holder Straight delicate jaws			5 1/8" (13 cm)
6006-04	SCANLAN® Dennis Delicate Needle Holder Straight jaws			7" (18 cm)
3003-422	SCANLAN® Microvascular Needle Holder Straight Jaws With Ratchet			7 1/8" (18.1 cm)



SCANLAN®

MICROSURGICAL NEEDLE HOLDERS

SCANLAN®

Micro Needle Holders

for 5-0 Suture



Cat. No.	Description	Handle	Recommended Suture Size	Total Length
3003-90A	SCANLAN® Titanium Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Streamline Box Lock with Ratchet	Round	5-0 and Smaller	4 3/4" (12 cm)
4155-06	SCANLAN® Webster Needle Holder Delicate DIAMOND DUST™ Jaws	Ring		5 1/8" (13 cm)
6006-24	SCANLAN® Micro Needle Holder DIAMOND DUST™ Jaws			
6006-50	SCANLAN® Castroviejo Micro Needle Holder Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	Round		5 3/4" (14.5 cm)
6006-56-3	SCANLAN® Castroviejo Micro Needle Holder Curved, 0.3mm tapered DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet			
3003-94-4	SCANLAN® Titanium Castroviejo Micro Needle Holder Curved, 0.4mm DIAMOND DUST™ Jaws Streamline Box Lock Without Ratchet			
6006-28	SCANLAN® Micro Needle Holder DIAMOND DUST™ Jaws	Ring	6 1/4" (16 cm)	



SCANLAN®
MICROSURGICAL FORCEPS

SCANLAN®
Super-Micro Fine Forceps



3303-410



3303-470

Cat. No.	Description	Handle	Total Length
NEW 3303-452	SCANLAN® <i>Super-Micro Fine</i> Titanium Forceps Angled 0.3 mm DIAMOND DUST™ Tips	Round	5" (12.7 cm)
NEW 3303-410	SCANLAN® <i>Super-Micro Fine</i> Titanium Forceps Straight 0.3 mm DIAMOND DUST™ Tips	Flat	
3303-470	SCANLAN® <i>Super-Micro Fine</i> Titanium Forceps Straight 0.3 mm DIAMOND DUST™ Tips	Round	7 1/4" (18.5 cm)
NEW 3303-472	SCANLAN® <i>Super-Micro Fine</i> Titanium Forceps Curved 0.3 mm DIAMOND DUST™ Tips		
NEW 3303-430	SCANLAN® <i>Super-Micro Fine</i> Titanium Forceps Straight 0.3 mm DIAMOND DUST™ Tips	Flat	7 1/4" (18.5 cm)
3303-434	SCANLAN® <i>Super-Micro Fine</i> Titanium Forceps Angled 0.3 mm DIAMOND DUST™ Tips	Flat	



SCANLAN®
MICROSURGICAL FORCEPS

SCANLAN®
Micro Forceps



3003-154-3

Cat. No.	Description	Metal	Handle	Total Length
4004-03-3A	SCANLAN® Micro Suture Tying Forceps 45 Degree Angled 0.3 mm, DIAMOND DUST™ Tips	Stainless Steel	Round	5 1/2" (14 cm)
3003-154-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Straight 0.3 mm DIAMOND DUST™ Tips	Titanium		6" (15 cm)
3003-156-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Curved 0.3 mm DIAMOND DUST™ Tips	Titanium		6" (15 cm)
3003-158-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Straight 0.3 mm DIAMOND DUST™ Tips	Titanium		7 1/4" (18.5 cm)
3003-159-3	SCANLAN® Micro Suture Tying Forceps Gerald Type Taper Curved 0.3 mm DIAMOND DUST™ Tips	Titanium		7 1/4" (18.5 cm)
4004-97-3A	SCANLAN® Micro Suture Tying Forceps 45 Degree Angled 0.3mm, DIAMOND DUST™ Tips	Stainless Steel		7 1/4" (18.5 cm)

SCANLAN®
Specialty Forceps



4004-29

Cat. No.	Description	Handle	Total Length
4004-05	SCANLAN® Jewelers Forceps Straight Sharp Smooth Tips	Flat	4 1/4" (11 cm)
4004-07	SCANLAN® Jewelers Forceps Straight Sharp DIAMOND DUST™ Tips		
4004-29	SCANLAN® Adson Forceps 1.2 mm Tips, 1x2 Teeth		4 3/4" (12 cm)
4004-19	SCANLAN® Jewelers Forceps Straight sharp DIAMOND DUST™ Tips		5 1/8" (13 cm)
4004-63	SCANLAN® McIndoe Forceps Serrated Tips, Fine 1x2 Teeth		6 1/4" (16.5 cm)



SCANLAN®

MICROSURGICAL SCISSORS

SCANLAN®

Super-Micro Fine Spring Style Scissors



Cat. No.	Description	Handle	Total Length
NEW 7007-712	SCANLAN® <i>Premier</i> ™ Super-Micro Fine Jacobson Scissors Curved Nano Blades	Flat	5" (12.7 cm)
NEW 7007-729	SCANLAN® <i>Premier</i> ™ Super-Micro Fine Jacobson Scissors Curved Nano Blades	Round	7 1/4" (18.5 cm)
NEW 7007-768	SCANLAN® <i>Premier</i> ™ Super-Micro Fine Jacobson Scissors Straight Nano Blades		
7007-124	SCANLAN® Jacobson Micro Scissors Curved Fine, Sharp Tips	Flat	7" (18 cm)
7007-120	SCANLAN® Jacobson Micro Scissors Straight Fine, Sharp Tips		
7007-140	SCANLAN® Jacobson Vannas Micro Scissors 45° Angled on Flat Micro Tips		



SCANLAN®

MICROSURGICAL SCISSORS

SCANLAN®

Micro Spring Style Scissors



Cat. No.	Description	Handle	Total Length
7007-404	SCANLAN® <i>Premier</i> ™ Castroviejo Micro Scissors 45° Angled Blades, 1 Blade with Bead	Flat	4 1/8" (10.5 cm)
7007-108	SCANLAN® Castroviejo Micro Scissors Curved Blades, Sharp Tips	Round	5 1/2" (14 cm)
7007-538	SCANLAN® <i>Premier</i> ™ Jacobson Micro Scissors Straight Sharp, Nano Blades	Flat	6" (15 cm)
7007-481	SCANLAN® Jacobson Micro Scissors 45° Angled Blades, 1 Blade with Bead tip	Flat	6 1/2" (16.5 cm)
7007-128	SCANLAN® Jacobson Micro Scissors Straight Blades, Sharp Tips	Round	7" (18 cm)
7007-130	SCANLAN® Jacobson Micro Scissors Curved Blades, Sharp Tips	Round	7" (18 cm)



SCANLAN®
MICROSURGICAL SCISSORS

SCANLAN®
Micro Ring Handle Scissors




Cat. No.	Description	Handle	Total Length
7007-282SC	SCANLAN® Iris-Eye Scissors Straight Super Cut™ Blades, Sharp Tips, Black Platinum Ring Handles	Ring	4" (10 cm)
7007-210SC	SCANLAN® Stevens Tenotomy Scissors Curved Super Cut™ Blades, Black Platinum Ring Handles		4 5/8" (11.8 cm)
7007-214SC	SCANLAN® Baby Metzenbaum Scissors Curved Super Cut™ Blades, Black Platinum Ring Handles		4 3/4" (12 cm)
7007-272SC	SCANLAN® Ragnell Dissecting Scissors Curved Super Cut™ Blades, Flat Tips, Black Platinum Ring Handle		6" (15 cm)
7007-292SC	SCANLAN® C.V.H. Scissors Straight Super Cut™ Blades, Sharp Tips, Black Platinum Ring Handles		6" (15 cm)
7007-211SC	SCANLAN® Stevens Tenotomy Scissors Curved Super Cut™ Blades Black Platinum Ring Handles		7" (18 cm)
7007-218-1SC	SCANLAN® Metzenbaum Scissors Curved Super Cut™ Blades, Fine Tips, Black Platinum Ring Handles		
7007-306SC	SCANLAN® Face-Lift Scissors Curved Super Cut™ Blades, Semi-Sharp Outer Blade, Black Platinum Ring Handles		7 3/4" (19.5 cm)



SCANLAN®
CLAMPS / HAEMOSTATS

SCANLAN®
CLAMPS / HAEMOSTATS




Cat. No.	Description	Total Length
9009-148	SCANLAN® Ochsner Forceps Ring Handle, 90° Angle, Very Fine Jaws	4 1/2" (11.5 cm)
5005-02	SCANLAN® Vascular (Mosquito) Forceps Ring Handle, Delicate Curved Jaws	4 3/4" (12 cm)
5005-64	SCANLAN® Bailey Forceps Ring Handle, 90° Angled Very Fine Jaws	7" (18 cm)



SCANLAN®
DILATORS

SCANLAN®
GRADUATED DILATORS



Cat. No.	Description	Total Length
2950-01	SCANLAN® Graduated Dilator Straight	5 1/2" (14 cm)
2950-02	SCANLAN® Graduated Dilator Curved	

SCANLAN®
MICRO VESSEL DILATORS



Cat. No.	Description	Total Length
4004-150	SCANLAN® Micro Vessel Dilator Spring Style, Flat Handle, Straight Tips	4 1/4" (11 cm)
4004-155	SCANLAN® Micro Vessel Dilator Spring Style, Round Handle, Straight Tips	4 7/8" (12.4 cm)



SCANLAN®
SINGLE-USE BULLDOG CLAMPS

SCANLAN®
VASCU-STATT® Bulldog Appliers



Cat. No.	Description
1001-506	SCANLAN® VASCU-STATT® Applier Straight, 8 1/2" (21.5 cm)
1001-506A	SCANLAN® VASCU-STATT® Applier Angled, 8 1/2" (21.5 cm)
1001-506M	SCANLAN® Micro VASCU-STATT® Applier Straight, 8 1/2" (21.5 cm)

SCANLAN®
VASCU-STATT® Bulldog Clamps



Cat. No.	Description	Clamping Pressure
1001-507	SCANLAN® VASCU-STATT® Approximator Mini, Straight (Red)	Regular
1001-508	SCANLAN® VASCU-STATT® Approximator Mini, Angled 45° (Red)	Regular
1001-540	SCANLAN® VASCU-STATT® II Approximator Straight (Green)	Firm
1001-499	SCANLAN® Micro VASCU-STATT® Micro, Straight (Purple)	Regular
1001-529	SCANLAN® Micro VASCU-STATT® II Micro, Straight (Black)	Firm



SCANLAN®

MICROSURGICAL Instrumentation

Performance, quality and exceptional tactile response during the most complex procedures.

For more information or to place an order,
please contact your local Scanlan representative,
or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

www.scanlaninternational.com

