


SCANLAN[®] *LEGACY*

High Precision Titanium
Instruments





Joseph
Scanlan

Russell
Scanlan

Tim
Scanlan

Adam
Scanlan

Chenoa
Scanlan

Dear Customer,

The Scanlan family has specialized in the design and manufacture of the highest quality surgical products since 1921. Now in our fourth generation, my family continues its dedication to providing you with the best surgical instrumentation available anywhere in the world.

Scanlan has continued this tradition with the development of the titanium LEGACY needle holders and forceps. These precision titanium instruments offer superior balance, handling and tactile feel. Now featuring the revolutionary SCAN-II® handle design that is unlike any other instrument available!

Sincerely,

Timothy M. Scanlan
President
Chief Executive Officer

Russel T. Scanlan
Chief Operating Officer


100 Years of Excellence... And Our Journey Continues

SCANLAN®

Features of LEGACY Titanium Forceps

Unique 10 mm Diameter and Extended Knurling

- Provides the titanium instruments with the optimized weight and balance.
- The extended knurling provides more surface area for the surgeon to grasp the instrument.



Perfect-Pyramid™
Provides a secure grasp.

Non-Magnetic

The properties of titanium prevent the needles from clinging to the instrument.



Never Shear™ Dual Guide™

Provides secure, reliable control without trauma to tissue or damage to needles.



SCAN-II® Handle

Revolutionary design that features a flat surface for control with round knurled sides for a secured grip and the ability to roll the instrument during procedures.



Round Handle

The original 10 mm round design provides the surgeon with a secure grip and the ability to roll the instrument during procedures.

Diamond-Dust™

Provides secure, reliable control without trauma to tissue or damage to needles.

Dual Guide™ Patent No. 8,147,512
The Never Shear™ Dual Guide™ design was developed in cooperation with John D. Puskas, MD

SCAN-II® Handle Design U.S. Patent No. D907,202 & D949,337





3003-942
Instrument is shown approximate size

SCANLAN®
LEGACY Micro Forceps

SCAN-II® Handle | *Never Shear™* Dual Guide™ | 10 mm, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No.	Description	Total Length
		3003-942	1.0 mm Straight DIAMOND DUST™ Platform Ring Tips	7 1/8" (18.1 cm)
		3003-944		8 1/4" (21 cm)
		3003-946		9" (23 cm)
		3003-948		10" (25.4 cm)



SCANLAN®
LEGACY Micro Forceps

Round Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No.	Description	Total Length
		3003-922	1.0 mm Straight DIAMOND DUST™ Platform Ring Tips	7 1/8" (18.1 cm)
		3003-924		8 1/4" (21 cm)
		3003-926		9" (23 cm)
		3003-928		10" (25.4 cm)



3003-974

SCANLAN®
LEGACY Micro Forceps

SCAN-II® Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No. 0.5 mm Tips	Cat. No. 1.0 mm Tips	Description	Total Length
		3003-962	3003-972	Straight DIAMOND DUST™ Platform Tips	7 1/8" (18.1 cm)
		3003-964	3003-974		8 1/4" (21 cm)
		3003-966	3003-976		9" (23 cm)
		3003-968	3003-978		10" (25.4 cm)



3003-914

SCANLAN®
LEGACY Micro Forceps

Round Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No. 0.5 mm Tips	Cat. No. 1.0 mm Tips	Description	Total Length
		3003-902	3003-912	Straight DIAMOND DUST™ Platform Tips	7 1/8" (18.1 cm)
		3003-904	3003-914		8 1/4" (21 cm)
		3003-906	3003-916		9" (23 cm)
		3003-908	3003-918		10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com





SCANLAN®
LEGACY Micro Forceps

SCAN-II® Handle | *Never Shear™* Dual Guide™ | 10 mm, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No.	Description	Total Length
		3003-952	1.0 mm Straight Tapered DeBakey Tips	7 1/8" (18.1 cm)
		3003-954		8 1/4" (21 cm)
		3003-956		9" (23 cm)
		3003-958		10" (25.4 cm)



SCANLAN®
LEGACY Micro Forceps

Round Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No.	Description	Total Length
		3003-932	1.0 mm Straight Tapered DeBakey Tips	7 1/8" (18.1 cm)
		3003-934		8 1/4" (21 cm)
		3003-936		9" (23 cm)
		3003-938		10" (25.4 cm)



SCANLAN®
LEGACY Micro Forceps

SCAN-II® Handle | *Never Shear™* Dual Guide™ | 10 mm, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No.	Description	Total Length
		3003-913	1.0 mm Straight Serrated Tips	7 1/8" (18.1 cm)
		3003-915		8 1/4" (21 cm)
		3003-917		9" (23 cm)
		3003-919		10" (25.4 cm)



SCANLAN®
LEGACY Micro Forceps

Round Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No.	Description	Total Length
		3003-903	1.0 mm Straight Serrated Tips	7 1/8" (18.1 cm)
		3003-905		8 1/4" (21 cm)
		3003-907		9" (23 cm)
		3003-909		10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com

3003-984



SCANLAN®
LEGACY Puskas™ Micro Forceps
SCAN-II® Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No. 0.5 mm Tips	Cat. No. 1.0 mm Tips	Description	Total Length
		3003-980		Straight DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	7 1/8" (18.1 cm)
		3003-984	3003-988		8 1/4" (21 cm)
3003-985			Curved DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	8 1/4" (21 cm)	
3003-986		3003-989	Angled DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	8 1/4" (21 cm)	
3003-990		3003-994	Straight DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	9" (23 cm)	
3003-992		3003-996		10" (25.4 cm)	

3003-778



SCANLAN®
LEGACY Puskas™ Micro Forceps
Round Handle | *Never Shear™* Dual Guide™ | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No. 0.5 mm Tips	Cat. No. 1.0 mm Tips	Description	Total Length
		3003-770		Straight DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	7 1/8" (18.1 cm)
		3003-774	3003-778		8 1/4" (21 cm)
		3003-775		Curved DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	8 1/4" (21 cm)
		3003-776	3003-779	Angled DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	8 1/4" (21 cm)
		3003-780	3003-784	Straight DIAMOND DUST™ Platform Tips 1 x 2 Micro Teeth	9" (23 cm)
		3003-782	3003-786		10" (25.4 cm)

3003-820



SCANLAN®
LEGACY Puskas™ Micro Forceps
SCAN-II® Handle | *Never Shear™* Dual Guide™ | Bayonet Style | 10 mm Diameter, Extended Knurling

Tip (enlarged to show detail)	Handle	Cat. No. 0.5 mm Tips	Cat. No. 1.0 mm Tips	Description	Total Length
	SCAN-II®	3003-820		Straight DIAMOND DUST™ Tips	8 3/4" (22 cm)

3003-810



SCANLAN®
LEGACY Micro Forceps
Round Handle | *Never Shear™* Dual Guide™ | Bayonet Style | 10 mm Diameter, Extended Knurling


Tip (enlarged to show detail)	Handle	Cat. No. 0.5 mm Tips	Cat. No. 1.0 mm Tips	Description	Total Length
	Round		3003-804	Straight DIAMOND DUST™ Platform Ring Tips	8 3/4" (22 cm)
		3003-816	3003-818	Straight DIAMOND DUST™ Platform Tips	
			3003-806	Straight Tapered DeBakey Tips	
		3003-810		Puskas™ Straight DIAMOND DUST™ Tips	

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com

Unique 10 mm Diameter and Extended Knurling

- Provides the titanium instruments with the optimized weight and balance.
- The extended knurling provides more surface area for the surgeon to grasp the instrument.



Perfect-Pyramid
Provides a secure grasp.

Non-Magnetic
The properties of titanium will prevent the needles from clinging to the instrument.

SCAN-II® Handle

Revolutionary design that features a flat surface for control with round knurled sides for a secured grip and the ability to roll the instrument during procedures.

Diamond-Dust™

Securely and delicately grips the needles to prevent slippage and any damage to the needle.

Round Handle

The original 10 mm round design provides the surgeon with a secure grip and the ability to roll the instrument during procedures.

3003-374



SCANLAN®
LEGACY Needle Holders
SCAN-II® Handle | 10 mm, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
		3003-372	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	7" (18 cm)
		3003-374			8" (20 cm)
		3003-376			9" (23 cm)
		3003-378			10" (25.4 cm)

3003-342

Instrument is shown approximate size



SCANLAN®
LEGACY Needle Holders
Round Handle | 10 mm Diameter, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
		3003-342	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	7" (18 cm)
		3003-344			8" (20 cm)
		3003-346			9" (23 cm)
		3003-348			10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com





SCANLAN®
LEGACY Needle Holders
SCAN-II® Handle | 10 mm, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
		3003-382	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	7" (18 cm)
		3003-384			8" (20 cm)
		3003-386			9" (23 cm)
		3003-388			10" (25.4 cm)



SCANLAN®
LEGACY Needle Holders
Round Handle | 10 mm Diameter, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
		3003-352	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	7" (18 cm)
		3003-354			8" (20 cm)
		3003-356			9" (23 cm)
		3003-358			10" (25.4 cm)



SCANLAN®
LEGACY Needle Holders
SCAN-II® Handle | 10 mm, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
		3003-332	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	3-0 and Smaller	7" (18 cm)
		3003-334			8" (20 cm)
		3003-336			9" (23 cm)
		3003-338			10" (25.4 cm)



SCANLAN®
LEGACY Jacobson Micro Needle Holders
Round Handle | 10 mm Diameter, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
		3003-362	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	3-0 and Smaller	7" (18 cm)
		3003-364			8" (20 cm)
		3003-366			9" (23 cm)
		3003-368			10" (25.4 cm)

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com

3003-381



SCANLAN®
LEGACY Puskas™ Jacobson Micro Needle Holders
Round Handle | 8 mm Diameter, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
	Round	3003-381	Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	7-0 and Smaller	8 1/4" (21 cm)
		3003-391			9 1/2" (24 cm)

3003-390



SCANLAN®
LEGACY Puskas™ Jacobson Micro Needle Holders
Round Handle | 8 mm Diameter, Extended Knurling

Jaws	Handle	Cat. No.	Description	Recommended Suture Size	Total Length
	Round	3003-390	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	5-0 and Smaller	9" (23 cm)

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com

Forceps



Platform Ring Tip	Total Length	Description	SCAN-II®	Round
	7-1/8" (18.1 cm)	1.0 mm Straight DIAMOND DUST™ Platform Ring Tips	3003-942	3003-922
	8-1/4" (21 cm)		3003-944	3003-924
	9" (23 cm)		3003-946	3003-926
	10" (25.4 cm)		3003-948	3003-928

Platform Tip	Total Length	Description	SCAN-II®	Round
	7-1/8" (18.1 cm)	0.5 mm Straight DIAMOND DUST™ Platform Tips	3003-962	3003-902
	8-1/4" (21 cm)		3003-964	3003-904
	9" (23 cm)		3003-966	3003-906
	10" (25.4 cm)		3003-968	3003-908
	7-1/8" (18.1 cm)	1.0 mm Straight DIAMOND DUST™ Platform Tips	3003-972	3003-912
	8-1/4" (21 cm)		3003-974	3003-914
	9" (23 cm)		3003-976	3003-916
	10" (25.4 cm)		3003-978	3003-918

DeBakey Tip	Total Length	Description	SCAN-II®	Round
	7-1/8" (18.1 cm)	1.0 mm Straight DeBakey Tips	3003-952	3003-932
	8-1/4" (21 cm)		3003-954	3003-934
	9" (23 cm)		3003-956	3003-936
	10" (25.4 cm)		3003-958	3003-938



Forceps



Serrated Tip	Total Length	Description	SCAN-II®	Round
	7-1/8" (18.1 cm)	1.0 mm Straight Serrated Tips	3003-913	3003-903
	8-1/4" (21 cm)		3003-915	3003-905
	9" (23 cm)		3003-917	3003-907
	10" (25.4 cm)		3003-919	3003-909

1x2 Teeth Platform Tip	Total Length	Description	SCAN-II®	Round
	7-1/8" (18.1 cm)	0.5 mm Straight 1X2 Teeth Platform Tips	3003-980	3003-770
	8-1/4" (21 cm)		3003-984	3003-774
	9" (23 cm)		3003-990	3003-780
	10" (25.4 cm)		3003-992	3003-782
	8-1/4" (21 cm)	1.0 mm Straight 1X2 Teeth Platform Tips	3003-988	3003-778
	9" (23 cm)		3003-994	3003-784
	10" (25.4 cm)		3003-996	3003-786
	8-1/4" (21 cm)	0.5 mm Curved 1X2 Teeth Platform Tips	3003-985	3003-775
	8-1/4" (21 cm)	0.5 mm Angled 1X2 Teeth Platform Tips	3003-986	3003-776
	8-1/4" (21 cm)	1.0 mm Angled 1X2 Teeth Platform Tips	3003-989	3003-779

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com



Bayonet Style Forceps



Platform Ring Tip	Total Length	Description	SCAN-II®	Round
1.0 mm	8-3/4" (22 cm)	1.0 mm Straight DIAMOND DUST™ Platform Ring Tips		3003-804

Platform Tip	Total Length	Description	SCAN-II®	Round
0.5 mm	8-3/4" (22 cm)	0.5 mm Straight DIAMOND DUST™ Platform Tips		3003-816
1.0 mm		1.0 mm Straight DIAMOND DUST™ Platform Tips		3003-818

DeBaKey Tip	Total Length	Description	SCAN-II®	Round
1.0 mm	8-3/4" (22 cm)	1.0 mm Straight DeBaKey Tips		3003-806

1x2 Teeth Platform Tip	Total Length	Description	SCAN-II®	Round
0.5 mm	8-3/4" (22 cm)	0.5 mm Straight 1X2 Teeth Platform Tips	3003-820	3003-810

Please check with your Scanlan Representative for Firm Sale status.
Custom modifications available upon request.

For a complete listing of models visit...
www.scanlaninternational.com



Needle Holders



7-0 Suture Size	Total Length	Description	SCAN-II®	Round
	7" (18 cm)	10 mm Knurling Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet	3003-372	3003-342
	8" (20 cm)		3003-374	3003-344
	9" (23 cm)		3003-376	3003-346
	10" (25.4 cm)		3003-378	3003-348
	8-1/4" (21 cm)	8 mm Knurling Straight DIAMOND DUST™ Jaws Regular Box Lock With Ratchet		3003-381
	9-1/2" (24 cm)			3003-391

5-0 Suture Size	Total Length	Description	SCAN-II®	Round
	7" (18 cm)	10 mm Knurling Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	3003-382	3003-352
	8" (20 cm)		3003-384	3003-354
	9" (23 cm)		3003-386	3003-356
	10" (25.4 cm)		3003-388	3003-358
	9" (23 cm)	8 mm Knurling Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet		3003-390

3-0 Suture Size	Total Length	Description	SCAN-II®	Round
	7" (18 cm)	Straight DIAMOND DUST™ Jaws Streamline Box Lock With Ratchet	3003-332	3003-362
	8" (20 cm)		3003-334	3003-364
	9" (23 cm)		3003-336	3003-366
	10" (25.4 cm)		3003-338	3003-368



LEGACY High Precision Titanium Instruments

For more information or to place an order,
 please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458
 International: 651-298-0997
 FAX: 651-298-0018
www.scanlaninternational.com