

A/V Access INSTRUMENTS



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921



A/V Access Instruments

Since 1921, Scanlan International, Inc. has been designing and manufacturing the highest quality surgical instrumentation. Our A/V Access instruments are a unique collection of instruments for procedures working within small to medium sized surgical fields. Each instrument is designed with the highest grade stainless steel or titanium alloys, and are handcrafted with meticulous precision... allowing excellent tactile response and control.


This collection has been selected from our full line of high quality surgical instrumentation (consisting of over 3000 designs). Please visit us at www.scanlaninternational.com to review additional selections and options.

TUNNELING SYSTEM

9009-14



9009-13



Cat. No.	Description	Length
9009-13	A/V Access Tuner Partial Curve	12" (30.5 cm)
9009-14	A/V Access Tuner	12" (30.5 cm)

9009-20



Interior Diameter



For grafts up to 8 mm

9009-22



Interior Diameter



For grafts up to 6 mm

Cat. No.	Description	Length
9009-20	A/V Access Tuner Sheath and Bullet Tip Large (Orange) for Grafts Up To 8 mm	12" (30.5 cm)
9009-22	A/V Access Tuner Sheath and Bullet Tip Small (Yellow) for Grafts Up To 6 mm	12" (30.5 cm)



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com



A/V Access Instruments

TUNNELING SYSTEM

9009-15



Cat. No.	Description	Length
9009-15	A/V Access Tunneler Full Curve	12" (30.5 cm)

9009-21



Interior Diameter



For grafts up to 8 mm

9009-23



Interior Diameter



For grafts up to 6 mm

Cat. No.	Description	Length
9009-21	A/V Access Tunneler Sheath and Bullet Tip Full Curve, Large (Orange) for Grafts Up To 8 mm	12" (30.5 cm)
9009-23	A/V Access Tunneler Sheath and Bullet Tip Full Curve, Small (Yellow) for Grafts Up To 6 mm	12" (30.5 cm)

9009-26



Cat. No.	Description	Length
9009-26	SCANLAN® V.I.P.™ Storage / Sterilization Tray	3.5" x 26" x 1.5" (8.9 cm x 66 cm x 3.8 cm)



A/V Access Instruments

FORCEPS

3003-922



Cat. No.	Description	Length
4004-76	Gerald Forceps, Straight Tips, Flat Handle 1 mm Serrated Tips	7" (18 cm)
4004-244	Gerald Forceps, Flat Handle 1 mm Straight DIAMOND DUST™ Ring Tips	7" (18 cm)
3003-912	Titanium SCANLAN® LEGACY Micro Forceps Never Shear™ Dual Guide™ 1.0 mm Straight Platform Tips, Round Handle	7 1/8" (18.1 cm)
3003-922	Titanium SCANLAN® LEGACY Micro Forceps Never Shear™ Dual Guide™ 1.0 mm Straight Platform Ring Tips, Round Handle	7 1/8" (18.1 cm)
3003-932	Titanium SCANLAN® LEGACY Micro Forceps Never Shear™ Dual Guide™ 1.0 mm Tapered DeBakey Tips, Round Handle	7 1/8" (18.1 cm)
4004-220	Dennis Micro Forceps, Ring Tips, Round Handle Straight DIAMOND DUST™ Tips, Counterbalanced	7 1/4" (18.5 cm)
4004-268	SCANLAN® DeBakey Needle Pulling Tissue Forceps, Round Handle Straight Micro Carbon-InLay Tips	7 1/4" (18.5 cm)
4004-43	DeBakey Forceps, Flat Handle 1 mm Straight Tips	7 3/4" (19.5 cm)

NEEDLE HOLDERS

3003-354



Cat. No.	Description	Length	Suture Size
6006-10	Fell Needle Holders 2 mm Straight Carbon-Inlay Jaws, Ring Handle (Gold Handle)	7" (18 cm)	3-0 and smaller
3003-462	Titanium SCANLAN® Micro Needle Holders Straight DIAMOND DUST™ Jaws, Ring Handle (Gold Handle)	7 1/4" (18.5 cm)	5-0 and smaller
6006-82	Jacobson Micro Needle Holders Straight DIAMOND DUST™ Jaws, Regular Box Lock, Round Handle	7 1/4" (18.5 cm)	5-0 and smaller
6006-170	Jacobson Micro Needle Holders Straight DIAMOND DUST™ Jaws, Regular Box Lock, Flat Knurled Handle	7 1/4" (18.5 cm)	5-0 and smaller
3003-354	Titanium SCANLAN® LEGACY Micro Needle Holders Straight DIAMOND DUST™ Jaws, Streamline Box Lock, Round Handle	8" (20 cm)	5-0 and smaller
3003-181	Jacobson Micro Needle Holders Straight DIAMOND DUST™ Jaws, Regular Box Lock, Round Handle	8" (20 cm)	7-0 and smaller

CLAMPS

5005-64



Cat. No.	Description	Length
5005-03	Vascular (Mosquito) Forceps Ring Handle, Curved Serrated Jaws	5" (12.7 cm)
2002-82	Dennis Multi-Purpose Clamp 4 cm 45° Angled Jaws	6 1/4" (16 cm)
5005-64	Bailey Forceps Ring Handle, 90° Angled Very Fine Serrated Jaws	7" (18 cm)



A/V Access Instruments

SCISSORS

7007-218SC



Cat. No.	Description	Length
7007-211SC	Stevens Tenotomy Scissors Curved Super Cut™ Blades Black Platinum Ring Handle	7" (18 cm)
7007-218SC	Metzenbaum Scissors Curved Super Cut™ Blades Black Platinum Ring Handle	7" (18 cm)
7007-252SC	Diethrich-Potts Scissors 45° Angled Super Cut™ Blades, Ring Handle	7" (18 cm)
7007-274SC	Ragnell Dissecting Scissors Curved Super Cut™ Blades Flat Tips, Black Platinum Ring Handle	7" (18 cm)
7007-442-8	Premier Jacobson Micro Scissors Curved 45° Angled Super Cut™ Blades, Flat Handle, Spring Style	8" (20 cm)

SUCTION

9009-960



Cat. No.	Description	Length
9009-960	Sundt™ Graduated Suction System , Rotatable, 3Fr., 4Fr., 5Fr., 6Fr., 7Fr., 9Fr., 11Fr.	5 1/4" (13.5 cm)

VASCULAR DILATORS / VESSEL PROBES

8008-460



Cat. No.	Description	Length
8008-456	SCANLAN® MEMORY Garrett Vascular Dilators / Vessel Probes 2 mm	5 1/2" (14 cm)
8008-458	SCANLAN® MEMORY Garrett Vascular Dilators / Vessel Probes 2.5 mm	5 1/2" (14 cm)
8008-460	SCANLAN® MEMORY Garrett Vascular Dilators / Vessel Probes 3 mm	5 1/2" (14 cm)
8008-462	SCANLAN® MEMORY Garrett Vascular Dilators / Vessel Probes 3.5 mm	5 1/2" (14 cm)
8008-464	SCANLAN® MEMORY Garrett Vascular Dilators / Vessel Probes 4 mm	5 1/2" (14 cm)

TEMPORARY OCCLUSION CLIPS (reusable)



Cat. No.	Description	Length
3795-01	Heifetz™ Temporary Occlusion Clip Appliers, Straight	7 1/4" (18.5 cm)
3795-14	Heifetz™ Temporary Occlusion Clip Straight, Firm Pressure, Gold 8 mm Jaw Length x 1.75 Jaw Width	
3795-18	Heifetz™ Temporary Occlusion Clip Curved, Firm Pressure, Gold 8 mm Jaw Length x 1.75 Jaw Width	
3795-20	Heifetz™ Temporary Occlusion Clip Straight, Light Pressure, Black 12 mm Jaw Length x 2.25 Jaw Width	
3786-23	Yasargil-type Temporary Occlusion Clip Angled 9 mm	
3786-24	Yasargil-type Temporary Occlusion Clip Slightly curved 9 mm	
3786-28	Yasargil-type Applying Forceps	7 1/8" (18.1 cm)




INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921



A/V Access Instruments

BULLDOG CLAMPS and APPLIERS



Cat. No.	Description
1001-506	VASCU-STATT® Applier Straight, 8 1/2" (21.5 cm)
1001-506M	Micro VASCU-STATT® Applier Straight, 8 1/2" (21.5 cm)
1001-499	Micro VASCU-STATT® Micro, Straight (Purple) Clamping Pressure 15-20 grams
1001-529	Micro VASCU-STATT® II Micro, Straight (Black) Clamping Pressure 30-40 grams
1001-550	VASCU-STATT® Plus Midi, Straight (White) Clamping Pressure 25-30 grams
1001-571	VASCU-STATT® II Plus Mini, Angled 45° (Green) Clamping Pressure 20-25 grams
1001-573	VASCU-STATT® II Plus Midi, Angled 45° (Yellow) Clamping Pressure 75-80 grams
1001-574	VASCU-STATT® II Plus Maxi, Straight (Orange) Clamping Pressure 165-175 grams


VEIN HOLDER

Cat. No.	Description
1001-760	SOLEM™ Vein Holder 6 1/2" (16.5 cm)





SILICONE LOOPS

Cat. No.	Description
1001-76	Red, Maxi
1001-77	Blue, Mini



TUNNELING SYSTEM



Cat. No.	Description	Length
2150-03	SCANLAN® V.I.P.™ Sterilization Tray Aluminum, Scan-Guard™ Silicone brackets & mat, Detachable Lid, Large capacity tray; Holds the entire SCANLAN® Tunneling System or 20 Instruments	29" x 8 3/8" x 1 1/2" (73.7 cm x 21.3 cm x 3.8 cm)

For more information or to place an order,
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458
International: 651-298-0997
FAX: 651-298-0018
www.scanlaninternational.com



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com

MK-B-avAccess Rev. B
© 2017 Scanlan International, Inc.