

SCANLAN® Minimally Invasive Dura Closure Set



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

SCANLAN®

Minimally Invasive Dura Closure Set

The SCANLAN® Minimally Invasive Dura Closure Set includes two extended needle holders (5 mm shaft), a knot pusher (5 mm shaft), a micro forceps and a tray. This allows for direct closure of the dura during minimally invasive spinal surgery, reducing the risk of cerebrospinal fluid leakage postoperatively.

Instrument Sterilization Tray

Cat. No.	Description
----------	-------------

2145-02	Dura Closure Set Metal Sterilization Tray, 10" x 16" x 2.5" (25.4 cm x 40.6 cm x 6.4 cm) Instrument capacity: 4-10 Firm Sale - May be available for evaluation
---------	---




SET NUMBER 9999-109

9009-360

Ring Handle Needle Holder


- Jaws open vertically
- DIAMOND DUST™ jaws
- Curved sideways shanks
- Working length 16.5 cm
- Ring handle
- 5 mm shaft
- 11 3/4" (30 cm)



9009-364

Jacobson Needle Holder

- Jaws open horizontally
- DIAMOND DUST™ jaws
- 30° angled down
- Working length 16.5 cm
- Round handle spring style
- 5 mm shaft
- 12 1/2" (32 cm)



9009-380


Chitwood Knot Pusher

- 30° angled down
- Working length 13.5 cm
- Round handle spring style
- 5.0 suture or smaller
- 5 mm shaft
- 11 3/4" (30 cm)



Additional Instruments

SUPER CUT™ Micro Scissors



- Straight or curved fine blades / blunt tips
- Working length 16.5 cm / 1 cm blade length
- 5 mm shaft diameter
- 10 1/2" (27 cm)

Cat. No.	Description
----------	-------------

9009-716SC	Straight fine blades, blunt tips, 1 cm blade length, total length, 10 1/2" (27 cm)
9009-718SC	Curved fine blades, blunt tips, 1 cm blade length, total length, 10 1/2" (27 cm)

Micro Forceps



- Round handle bayonet style
- Straight 0.8 mm DIAMOND DUST™ tips
- 10" (25.4 cm)

Cat. No.	Description
----------	-------------

3789-20	0.8 mm tips, 10" (25.4 cm)
---------	----------------------------

For more information or to place an order, please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

www.scanlaninternational.com



SCANLAN ★

INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlaninternational.com

MK-B-DURA Rev. F
© 2017 Scanlan International, Inc.