



NEW Titanium Dennis Multi-Purpose Clamps



Cat. No.	Description	Length
3003-280	Titanium, straight 2 x 3 DeBakey jaws	8 1/2" (21.5 cm)



Cat. No.	Description	Length
3003-282	Titanium, 45° angled 2 x 3 DeBakey jaws	8" (20 cm)




Cat. No.	Description	Length
3003-286	Titanium, 90° angled 2 x 3 DeBakey jaws	5 1/2" (14 cm)




INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921


NEW Stainless Steel Dennis Multi-Purpose Clamps



Cat. No.	Description	Length
2002-280	Stainless steel, straight 2 x 3 DeBakey jaws	8 1/2" (21.5 cm)



Cat. No.	Description	Length
2002-282	Stainless steel, 45° angled 2 x 3 DeBakey jaws	8" (20 cm)



Cat. No.	Description	Length
2002-286	Stainless steel, 90° angled 2 x 3 DeBakey jaws	5 1/2" (14 cm)

For more information or to place an order, please contact your local Scanlan representative, or call...

U.S. & Canada:
800-328-9458

International:
651-298-0997

FAX: 651-298-0018
www.scanlaninternational.com



INTERNATIONAL