

Extracranial-Intracranial (EC/IC) Bypass Instrumentation



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921



Committed To Our Customers

Dear Customer,

The Scanlan family has specialized in the design and manufacture of the world's finest surgical instrumentation and single-use products since 1921. Now in our fourth generation, our family continues its dedication to providing you with the best surgical products available in the world. This brochure is a collection of our precision products, all designed specifically to address your needs within the operating room.

Our goal is to continuously enhance our designs to serve you, your hospital and ultimately, your patients. Thank you for your continued support and demand for the finest surgical products.

Sincerely,

Timothy M. Scanlan
President and Chief Executive Officer

Russell T. Scanlan
Chief Operating Officer



Adam
Scanlan

Tim
Scanlan

Russell
Scanlan


Chenoa
Scanlan

EC / IC Bypass Instrumentation Set

Three lengths of instruments have been developed to accommodate the requirements for the varying depths of cranial bypass.

Additional features of this offering include:

- A precise 0.3 mm tip
- Scanlan's exclusive Diamond Dust™ to control needle rotation and enable atraumatic handling of needles, suture and tissue
- Streamlined tips on all forceps and needle holders provide enhanced visibility of micro needle and sutures
- Titanium forceps eliminate needle magnetism
- Metal instrument trays prevent damage to instruments during transport, cleaning and sterilization



0.3 mm 45° Tying Forceps



0.3 mm Gerald Type Taper
Straight Forceps



0.3 mm Gerald Type Taper
Curved Forceps



1 mm Ring Tip
Forceps



0.3 mm Curved
Needle Holder



45° Vannas
Micro Scissors



SCANLAN®

EC / IC Bypass Instrumentation Set

Short Instrument Set 9999-162

SCANLAN®
Forceps
Round Handle

Instrument is shown approximate size



SCANLAN® Micro Suture Tying Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-154-3	0.3 mm Straight Platform Gerald Type Taper DIAMOND DUST™	15.0 cm (6.00")

Instrument is shown approximate size



SCANLAN® Micro Suture Tying Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-156-3	0.3 mm Curved Platform Gerald Type Taper DIAMOND DUST™	15.0 cm (6.00")

Instrument is shown approximate size



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-164	1 mm Straight DIAMOND DUST™ Ring Tip	15.0 cm (6.00")

Instrument is shown approximate size



SCANLAN® Micro Suture Tying Forceps, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	4004-03-3A	0.3 mm Angled 45° Angle DIAMOND DUST™ Tips	14.0 cm (5.50")



SCANLAN®

EC / IC Bypass Instrumentation Set

Short Instrument Set 9999-162

SCANLAN®

Needle Holders

Round Handle

Instrument is shown approximate size



SCANLAN® Castroviejo Micro Needle Holder, Stainless Steel

Qty.	Cat. No.	Description	Recommended Suture Size	Total Length
1	6006-56-3	0.3 mm Tapered Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	14.5 cm (5.75")

SCANLAN®

Jacobson Vannas Micro Scissors

Flat Handle

Instrument is shown approximate size



SCANLAN® Jacobson Vannas Micro Scissors, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	7007-140	45° Angle Micro Fine Blades	18.0 cm (7.00")

See pages 10-11 for optional SCANLAN® V.I.P.™ Sterilization Trays



SCANLAN®

EC / IC Bypass Instrumentation Set

Medium Instrument Set 9999-163

SCANLAN®
Forceps
Round Handle

Instrument is shown approximate size



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-158-3	0.3 mm Straight DIAMOND DUST™ Gerald Type Taper Platform Tips	18.5 cm (7.25")

Instrument is shown approximate size



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-159-3	0.3 mm Curved DIAMOND DUST™ Gerald Type Taper Platform Tips	18.5 cm (7.25")

Instrument is shown approximate size



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-164	1 mm Straight DIAMOND DUST™ Ring Tips	15.0 cm (6.00")

Instrument is shown approximate size



SCANLAN® Micro Forceps, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	4004-97-3A	0.3 mm Angled DIAMOND DUST™ 45° Angle	14.0 cm (5.50")



SCANLAN®

EC / IC Bypass Instrumentation Set

Medium Instrument Set 9999-163

SCANLAN®

Needle Holders

Round Handle

Instrument is shown approximate size



SCANLAN® Jacobson Micro Needle Holder, Titanium

Qty.	Cat. No.	Description	Recommended Suture Size	Total Length
1	3003-177-3	0.3 mm Tapered Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	18.5 cm (7.25")

SCANLAN®

Jacobson Vannas Micro Scissors

Flat Handle

Instrument is shown approximate size



SCANLAN® Jacobson Vannas Micro Scissors, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	7007-140	45° Angle Micro Fine Blades	18.0 cm (7.00")

See pages 10-11 for optional SCANLAN® V.I.P.™ Sterilization Trays



SCANLAN®

EC / IC Bypass Instrumentation Set

Long Instrument Set 9999-164

SCANLAN®

Forceps

Bayonet Style | Flat Handle



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-842	0.3 mm Straight DIAMOND DUST™	22.0 cm (8.75")



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-841	0.3 mm 45° Angle DIAMOND DUST™	22.0 cm (8.75")



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-840	0.3 mm Curved Up DIAMOND DUST™	22.0 cm (8.75")



SCANLAN® Micro Forceps, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	3003-851	1 mm Straight DIAMOND DUST™ Ring Tip	22.0 cm (8.75")



SCANLAN®

EC / IC Bypass Instrumentation Set

Long Instrument Set 9999-164

SCANLAN®

Micro Needle Holders

Bayonet Style | Round Handle



SCANLAN® Micro Needle Holder, Stainless Steel

Qty.	Cat. No.	Description	Recommended Suture Size	Total Length
1	3789-14-3	0.3 mm Curved Up DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	23.0 cm (9.00")

SCANLAN®

Micro Scissors

Bayonet Style | Flat Handle



SCANLAN® Micro Scissors, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	3788-02	Vannas 45° Angled Up On Flat Micro Fine Blades	22.0 cm (8.75")

See pages 10-11 for optional SCANLAN® V.I.P.™ Sterilization Trays



Optional

SCANLAN® V.I.P.™ Sterilization Trays



2138-03

Description

Silver-gray color | Scan-Guard™ Silicone Inserts
25.4 x 25.4 x 3.8 cm (10.00" x 10.00" x 1.50")
Capacity: 8-10

2140-03

Description

Silver-gray color | Scan-Guard™ Silicone Inserts
25.4 x 38.0 x 3.8 cm (10.00" x 15.00" x 1.50")
Capacity: 10-12



Optional
SCANLAN® V.I.P.™ Sterilization Trays



2141-03

Description

Silver-gray color | base with accessory compartment | Scan-Guard™ Silicone Inserts
25.4 x 38.0 x 3.8 cm (10.00" x 15.00" x 2.6")
Capacity: 15+



EC / IC Bypass Instrumentation Set

For more information or to place an order,
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

www.scanlaninternational.com



SCANLAN ★

Surgical products designed and manufactured by the Scanlan family since 1921.
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlangroup.com

USA



MK-B-ECIC Rev. C